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
Fax. ++ 49 (0) 6261 – 9770-69

GENERAL ITEMS

These operating instructions provide the necessary information required to operate and maintain the sealer.

It is in your interest to make your operating personnel aware of these instructions prior to using the device.

The HAWO packaging machines are welding and sealing devices for the welding of thermoplastic films that means duplexing films according to the heatwelding and heatsealing method.

The special characteristics are described more detailed in „appropriate use of the device“.

These sealing machines are intended for industrial and business use and should only be used for the recommended purpose.

The equipment should only be used by a skilled person. The machines are only intended for use in dry surroundings and must be protected from humidity.

Before cleaning: Disconnect from the mains! Clean the machine only with a dry or damp soft cloth and a mild cleaning agent. Do not allow any water to find its way into the machine. Caution!

Never wet clean the machine!

In the interest of technological progress, we reserve the right to modify the machine.

As all other technical devices, also our devices are wearing out. In order to warranty a faultless functioning, the devices have to be opened by an experienced person at least once a year and have to be cleaned thoroughly and checked re-Safety.

The machine is only to be opened after the power plug has been pulled out.

Please observe these instructions very carefully.

If the machine is used by a person who is not able to read the operation instructions, the operator is obliged to undergo training.

Should the customer make use of the warranty, we reserve the right to make technical modifications in the interest of technical progress.

Regarding spare parts and deterioration parts please only use parts recommended and approved by HAWO.

Hawo will not accept any responsibility for resulting damages.

QUALITY TESTING

Every machine is subject to constant quality testing of the individual parts and of the finished assembly.

We demand this quality testing from our suppliers, too. We cannot, however, warranty the careful execution of their testing.

All machines are built in accordance with the relevant regulations.

Each electrical device is controlled according to the VDE-regulations.

DELIVERY AND CARRIAGE

Every machine is supplied in its original packaging.

Please check the contents carefully.

Before disposing the packaging material, please set machine into operation to make sure that no parts are left in the packaging.

Action to be taken when damaged during transit (even when the packaging is undamaged):
Damage which occurred during transit and which is not discovered under normal circumstances or which is reported too late, cannot be recovered.

Therefore please check the consignment immediately!

Should the consignment be damaged on the outside, a complaint should be made on acceptance and the damage reported (copy to supplier).

If damage is discovered after unpacking, a statement of the facts or negotiated letter should be sent immediately to the carriers. In this case, the packaging must be kept.

The following periods for notifying damage should be observed for ensuring a claim against the carrier:

Railway	:	within 7 days
Forwarding agent	:	within 4 days
Mail	:	within 24 hrs

In the event of damage an acknowledgement from the carrier on type and cause of damage is necessary. All other damage, for which one of the above mentioned occurrences cannot be proved, is attributable to the consignee or goes against transport insurance, if such was taken.

INSTALLATION

The device can be installed in each normal work-room.

Please be careful to keep the machines away from dust, steam and splashed water. These can cause the machine to malfunction.

The device must not be assembled and used in explosion-hazardous areas.

ELECTRICAL CONNECTION

CAUTION!

The machine is only to be connected after you are sure of the conformity of the mains voltage and the voltage figures on the machine type plate.

Only the switched-off machine may be connected with the supply mains.

CAUTION!

In case of any repair on the machine or on the power supply unit, first of all, pull out the power plug. Please only use plug connectors and plugs that are properly grounded.

WARRANTY

We grant the legal warranty for our sealers, based on an eight-hour duty. This means, in case of eventual material or production deficit we either deliver free replacement or repair the concerned parts free of charge (provided that you returned it carriage and we checked the machine thoroughly).

Regarding outside products and parts we can only transmit the warranty granted to us by the manufacturer of these parts. All those parts liable to wear and tear, like heater, PTFE-ribbon, rubber and plastic parts, such as seals, glass and relays are excluded from warranty.

Costs resulting from fitters travel are excluded from warranty.

Further claims or rights of recourse of third parties are not acknowledged.

The total warranty period is not being extended after warranty servicing has been done.

A case of warranty does not justify the non-payment, even partly, of our invoices.

Work only can be done under warranty if the following conditions are fulfilled:

The machine was used in workmanlike fashion and in accordance with the operating instructions.

Complaints were notified immediately, at the very latest 8 days after notice.

No actions or modifications without our permission have been carried out at the machine. The machine must be returned free of charge to the manufacturer.

The customer must have fulfilled his obligations.

Transport damages are excluded from warranty.

If the machine is damaged on return as a result of inadequate packaging the warranty claim lapses.

Necessary repairs of damages caused by an improper use of the machine and no more under warranty, will be effected immediately at our best prices.

In the interest of technological progress, we reserve the right to modify the machine.

We, hawo GmbH, Industriegebiet, D-74847 Obrigheim, do not accept any responsibility for resulting damages.

SAFETY OF DEVICES

The device has the GS test mark.

GENERAL SPECIFICATIONS

The welding device hpl ISZ is a portable small device designed for thermoplastic film via the heat impulse method.

The welding device consists of the impulse power supply unit "WING", the welding tong with a special cable and plug for the connection with the power supply unit.

The machine is only meant for commercial and industrial use and may only be used for the stipulated intended use with thermoplastic film.

The machine is not suitable for packing food and hazardous material. These may only be used after prior approval of the manufacturer.

One person is necessary to operate the machine.

INFORMATION REGARDING NOISE INTENSITY

The device works very silently.

INSTALLATION

The welding unit hpl ISZ is intended only for use in normal and dry surroundings.

SETTING INTO OPERATION

1. Put mains connection plug of power supply into socket (pay attention to the voltage indication on machine type plate)
2. Connect welding tong with socket power supply unit by plugging in the special plug.
3. Turn control knob "Welding Time" at power supply unit to position "0"
4. Switch on power supply unit by pressing tumbler to position "0"
5. Adjust welding time by turning the control knob "WeldingTime" to the right.

CAUTION!

The welding power may only be regulated step by step from 0 to 10. It depends on the kind of film which is to be welded and on the thickness of the film, therefore trial welding is necessary.

6. Insert material between the welding beams and close the tong by manual operation.

CAUTION!

Do not put any hard objects between the welding beams, because they may damage the heating ribbon as well as the PTFE glass tissue; the density of the welding seam in this spot could not be guaranteed anymore.

As soon as the welding beams are closed, the welding impulse is released and the welding process is performed according to the adjusted time value. The welding period is indicated by a buzzer signal in the device.

7. After the rundown of the welding impulse (buzzer signal ceases) the tong must be closed for another 2-5 sec. to secure the necessary re-cooling respectively the stability of the seam.
8. Open tong and take the films out.

9. Check seal quality :
 - A. Welding is too weak, seam can be opened easily- welding time is too short
 - B. Welding is too strong, seam is scorched – welding time is too long
10. After finishing welding operation switch off power supply unit by pressing tumbler switch into position "O".

REPLACEMENT OF DETERIORATION PARTS

CAUTION!

These jobs should only be carried out by skilled electricians. In case of any maintenance on the machine first of all pull the mains connection plug.

To keep up the faultless function of your device and to avoid unnecessary damages, we recommend the exclusive use of hawo- genuine parts.

PTFE-ribbon, heating ribbon and silicon rubber:

1. Dismount welding beam, remove stamp.
2. Pull off the defective PTFE-ribbon.
3. Loosen the M4 x 6- screws at both ends of the welding beam and remove the heating ribbon.
4. Replace silicon rubber, if necessary.
5. Insert new heating ribbon. Both ends of the heating ribbon must be double. The end shows to the ram body.
6. Fix heating ribbon by the M4 x 6- screws.
7. Cut new PTFE-ribbon true to size (if using meter storage), pull off the protection film and stick the self-adhesive PTFE-ribbon onto the welding ram.
8. Remount welding ram; pay attention to parallelity.
9. Fix again connection cable.

MAINTENANCE

As all other technical devices, also our devices are wearing out. In order to guarantee a faultless functioning, the devices have to be opened by an experienced person at least once a year and have to be cleaned thoroughly and checked re-safety.

We recommend to control the holding support of the operation lever every three months and to put some lubricating oil.

In case of spare- or deterioration parts orders please state the following data

- | | | |
|----|----------------------|---|
| 1. | Machine type | on machine type plate |
| 2. | Machine number | on machine type plate |
| 3. | Year of construction | on machine type plate |
| 4. | Spare part number | according to included
spare parts list |

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hawa GmbH
hawa GmbH, Obere Au 2, D-74847 Obrigheim, Germany