de	<b>Gebrauchsanweisung</b> Gasinsufflationssystem
en	<b>Instructions for use</b> Gas insufflation system
fr	<b>Mode d'emploi</b> Système d'insufflation de gaz
it	<b>Istruzioni per l'uso</b> Sistema di insufflazione di gas
es	<b>Instrucciones de uso</b> Sistema de insuflación de gas
pt	<b>Instruções de utilização</b> Sistema de insuflação de gás
cs	<b>Návod k použití</b> Systém pro insuflaci plynu
nl	<b>Gebruiksaanwijzing</b> Gasinsufflatiesysteem
pl	<b>Instrukcja użycia</b> System insuflacji gazu

<input type="checkbox"/> REF	05.000.22.568	<input type="checkbox"/> REF	05.000.22.581
<input type="checkbox"/> REF	05.000.22.574	<input type="checkbox"/> REF	05.000.22.582
<input type="checkbox"/> REF	05.000.22.575	<input type="checkbox"/> REF	05.000.22.583
<input type="checkbox"/> REF	05.000.22.580	<input type="checkbox"/> REF	05.000.22.584

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## Gas insufflation system

### Intended use

The gas insufflation system is used for applying and filtering gas (CO<sub>2</sub>) during laparoscopic and endoscopic procedures, e.g. inducing a pneumoperitoneum using an insufflation pump (insufflator) on a suitable patient connection such as a Ver-ess needle or a trocar.

### Intended users

The device may only be used by trained medical personnel.

### Patient group

The device is intended for use on adults and children.

### Description

This product is an accessory for a medical device. (CO<sub>2</sub> insufflator).

The system is made up of connecting tubes with appropriate connections for the particular insufflation equipment and filter. The filter protects the patient against any contamination in the insufflation gas and the insufflation pump against back-flow from the patient.

### Technical performance characteristics

REF	Hose length [cm]	Filter	Patient connection	Device connection
05.000.22.568	300 and 15	yes	Luer lock (EN ISO 80369-7)	Shortened funnel
05.000.22.574	300	yes	Luer lock (EN ISO 80369-7)	Filter 15M / 22F (EN ISO 5356-1)
05.000.22.575	300	yes	Luer lock (EN ISO 80369-7)	Filter 15M / 22F (EN ISO 5356-1)
05.000.22.580	300	yes	Luer lock (EN ISO 80369-7)	Filter with stepped connector
05.000.22.581	300 and 30	yes	Luer lock (EN ISO 80369-7)	Smooth end
05.000.22.582	300 and 15	yes	Luer lock (EN ISO 80369-7)	Shortened funnel
05.000.22.583	300 and 10	yes	Luer lock (EN ISO 80369-7)	Filter 15M (EN ISO 5356-1)
05.000.22.584	300 and 15	yes	Luer lock (EN ISO 80369-7)	Luer lock (EN ISO 80369-7)

### Contraindications

There are no known contraindications for using the gas insufflation system.

### Combination

The gas insufflation system has standard male Luer lock connectors in accordance with EN ISO 80369-7 on the patient end of the hose, which can be used to form a secure connection with products such as Veress needles or with female Luer lock connectors.

Depending on the product (REF), there are different connectors on the device connection end of the gas insufflation system to ensure compatibility with as many devices as possible. A suitable gas insufflation system must be chosen depending on the insufflation pump connection (see section "Performance characteristics" Device connection).

### Direction for use

A basic understanding of the insufflation technique as well as laparoscopic/endoscopic surgery is a prerequisite for using this system.

! Open the packaging under aseptic conditions and transfer to the instrument nurse in the sterile field.  
! Transfer the machine tube end (with filter) from the sterile field to the non-sterile circulating nurse and connect to the insufflation pump.

! Connect the patient tube end to the Luer Lock connection on the insufflation needle or the trocar.  
! Ensure that the connections and joints are tight to prevent loss of gas or pressure. The gas insufflation system must not be kinked because this can lead to incorrect displays (flow and pressure) on the insufflation pump.

! Start the gas insufflation on the insufflation pump according to the instructions given by the surgeon.

### Duration of application

The gas insufflation system can be used for short term application over the period of the surgery.

### Safety information

This product is intended for single use only. Sterilisation and cleaning (for re-use) can impair the function and mode of operation of the product and lead to a hazard for the patient.

When using the product, the clinical risks posed by the CO<sub>2</sub> gas must be observed:

! Hypothermia

! Acidosis as a result of hypercapnia

Ensure that the packaging is undamaged before using the product. The product must not be used under any circumstances if the sterile packaging has been damaged or unintentionally opened, as sterility can no longer be guaranteed. Check the use-by date before use. The product must not be used after the use-by date has expired.

The gas insufflation system must be stored in a dry place in their original packaging at room temperature and out of direct sunlight. If the storage conditions are not observed, the product properties may be affected. Should the packaged products have been exposed to environmental conditions exceeding those specified for a longer period of time, then the products may no longer be used.

### Disposal

Dispose of the product in accordance with your internal instructions and local regulations.

### Additional information

All serious incidents that occur in connection with the product must be reported to the manufacturer and the competent authorities.

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P. J. Dahlhausen & Co. GmbH
Dahlhausen Medizintechnik GmbH