
TUBE HOLDER

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Pictogram Legend

Pictogrammes listed below you'll find on the product packaging, if applicable.

Specially for children

MRT suitable

Catalogue number

Batch code

Use by

Content (in pieces)

Consult instructions for use

Keep away from sunlight

Keep dry

Do not reuse

Single patient use

Medical devices

CE marking

Manufacturer

TUBEHOLDER

en

FOREWORD

These instructions for use are valid for tube holders.

These instructions for use serve to inform the physician, nursing staff and patient/user in order to ensure technically correct handling.

When used for the first time, the products must be selected, used and inserted by appropriately trained physicians or medical professionals/medical device consultants. Please read the instructions for use carefully before using the product for the first time! Please keep the package for as long as you use the product.

It contains important information on the product!

1. INTENDED PURPOSE

Tube holders are used to fix and secure tracheostomy tubes in the tracheostoma of tracheotomised or laryngectomised patients.

2. SAFETY INFORMATION

The tube holder is a single-patient product and is thus intended exclusively for use - including repeated use - in one and the same patient.

A secure hold must always be monitored to prevent accidental loosening.

To prevent infections, replace soiled or wet tracheostomy tube straps. See hygiene instructions.

The tracheostomy tube strap must be selected and adjusted according to the neck width.

If too loose, this can lead to injury to the tracheostoma and dislocation of the tracheostomy tube.

If the tracheostomy tube is set too tight, pressure injuries may occur.

The tubeholder should not be fitted too tightly around the neck.

There should ideally be a distance of approx. 1 cm between tube holder and neck.

To prevent the tracheostomy tube from coming off by accident, the tube holder must be checked regularly to make sure it is still securely fixed in place.

The following applies for use with children: The tracheostomy tube strap must only be used in a supervised environment.

A secure hold must always be monitored to prevent accidental loosening.

3. INDICATION

The product is indicated for securing the tracheostomy tube in tracheotomized or laryngectomized patients.

4. CONTRAINDICATION

Do not use if the patient is allergic to the material.

5. PRODUCT DESCRIPTION

The tube holders feature various locking mechanisms.

- “H” variants: with plastic hook
- “K” variants: with velcro fastener
- “HM” variants: with metal hook

Tube holders are available in white, beige, blue and transparent as well as coloured (for children).

The tube holders consist of a foam strap surrounded by soft, non-woven fabric.

KACLIP® beige and **KACLIP®** white are made of a soft, elastic textile material, **KACLIP®** Clear H is made of an elastic polyurethane.

Tube holders are available in different lengths and can be adjusted to different neck sizes.

Please see the size table in the appendix.

6. APPLICATION

Attach the two parts of the tube holder to each other by means of the adherent fastener, so that the correct length (neck size) is obtained.

WARNING!

The tracheostomy tube strap must be selected and adjusted according to the neck width.

If too loose, this can lead to injury to the tracheostoma and dislocation of the tracheostomy tube.

If the tracheostomy tube is set too tight, pressure injuries may occur.

Place the tube holder around your neck with the hooks/hook-and-loop fasteners facing outward.

Attach the two hooks/hook-and-loop fasteners to the lateral eyelets on both sides of the neck flange.

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OPTIFLAUSCH®:

Adjust the short fastener tape with the stretch insert to the required neck size and fasten it at the fleecy side of the holder.

If the **OPTIFLAUSCH®** the holder is too long it can be shortened with a pair of scissors to the required length.

7. HYGIENE INSTRUCTIONS

To prevent infections, replace soiled or wet tracheostomy tube straps.

NECKFIX®:

Cleaning, disinfection, re-sterilisation and re-use can impair the safety and function of the product and are therefore not permissible!

KACLIP®, OPTIFLAUSCH®:

The products must not be cleaned or disinfected.

KACLIP® CLEAR:

The product may be cleaned by hand washing up to **40° C**.

Do not wring or spin-dry.

The tube holders must be completely dry before reuse.

Do not dry in a dryer.

SUPRAFIX®:

Hooks and hook-and-loop fasteners must be fastened prior to washing.

The tube holders may only be washed by hand up to **30 °C**.

Do not wring or spin-dry.

The tube holders must be completely dry before reuse.

Do not dry in a dryer.

8. STORAGE

Store in a dry environment away from sunlight and/or heat.

9. SERVICE LIFE

NECKFIX®:

The maximum period of use / wearing period is **24** hours.

KACLIP® und OPTIFLAUSCH®:

The tube holder is a single-patient product and also suitable for repeated use.

The service life depends on the individual frequency of use as well as the condition and on individual handling of the product.

When changing tubes, the tube holder in use must be placed temporarily on a clean surface.

Before being reattached, the product must be checked for visible damage and for functionality.

The product may not be reused if it is damaged, heavily soiled or heavily soaked, as this could jeopardise safe use.

In this case, the product must be replaced immediately.

The tube holder must be disposed of and replaced with a new tube holder no later than 7 days after first being used.

SUPRAFX®:

The tube holders are intended for repeated use and can be cleaned (washed).

The service life depends on the individual frequency of use and on individual handling

After 7 days at the latest, the tube holders must be disposed of.

The products must be replaced immediately when damaged.

10. DISPOSAL

The product must only be disposed of in accordance with the national regulations.

11. LEGAL NOTICES

The manufacturer Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for functional deficiencies, injuries, infections, and/or other complications or adverse events caused by unauthorised product alterations or improper use, care, and/or handling.

In particular, Andreas Fahl Medizintechnik-Vertrieb GmbH will not accept any liability for damages caused by modifications to the product or by repairs, if these modifications or repairs were not carried out by the manufacturer himself.

This applies – to the extent permitted by law – both to damages to the products themselves caused thereby as well as to any consequential damages caused thereby.

When using the product beyond the period of use specified under Section 9, and/or use, treatment, maintenance (cleaning, disinfection) or storage of the product in non-compliance with the specifications laid down in these instructions for use, Andreas Fahl Medizintechnik-Vertrieb GmbH will, as far as legally permissible, be free of any liability, including liability for defects.

Should a serious incident occur in connection with this product of Andreas Fahl Medizintechnik-Vertrieb GmbH, then this must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient are domiciled.

Sale and delivery of all Andreas Fahl Medizintechnik-Vertrieb GmbH products is carried out exclusively in accordance with our General Terms and Conditions of Business (GTC) which can be obtained directly from Andreas Fahl Medizintechnik-Vertrieb GmbH.

The manufacturer reserves the right to make changes to the product at any time without notice.

KACLIP®, SUPRAFX®, NECKFIX® und OPTIFLAUSCH® are registered trademarks of Andreas Fahl Medizintechnik-Vertrieb GmbH, Cologne in Germany and the EU member states.

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