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

 **Read the instructions for use carefully before application!****Product description:**

The PNEUMOCATH® Drainage Catheters consist of a radiopaque PE tube with ring markings and side openings for suction, inserted into a puncture cannula with attached sterile sheath. At the proximal end, a three-way-stopcock with step adapter is firmly connected by means of a metal thorn. The catheters are inserted through the needle. Due to its small size and easy handling, the PNEUMOCATH® thoracic drainage is also suitable for paediatrics. The NEO-PNEUMOCATH® models offer higher flow rates and better suction results for drainage due to larger lumens. The Drainage Kits can contain various components that are considered necessary for the surgically invasive procedure. These models are also commonly grouped and referred to as PNEUMOCATH®.

**Packaging content:**

The contents of the packaging and the specifications of the individual components are indicated on the packaging and the label.

**Valid for:**

All PNEUMOCATH® models (Drainage Catheters and Drainage Kits).

**Indications and Medical Purpose:**

The PNEUMOCATH® Drainage Catheters and Drainage Kits are used for suctioning unwanted effusions of all kinds from the pleural cavity (pleural effusion (e.g. haemothorax, pleural empyema, chylothorax)), for removing air or gases from the pleural cavity (pneumothorax (spontaneous, traumatic, iatrogenic)) or for pleurodesis. Pleural puncture (thoracocentesis) can be used for diagnostic or therapeutic purposes. The aims of the application are to initiate further therapeutic measures and, among other things, significantly to restore the physiological pressure conditions of the pleural cavity and to relieve symptoms or, in the case of pleurodesis, to obliterate the pleural leaves in order to prevent a recurrence of effusion.

 **Application duration:**

The product is only suitable for short-term (< 30 days) use.

 **Contraindications:**

There are no absolute contraindications, especially in case of vital (life-threatening) indication. It is at the discretion of the attending Healthcare Professional to assess this.

 **Relative contraindications:**

*It is necessary to weigh the benefit against the risk by the user!*

Unwilling uncooperative patient; Coagulopathy / haemorrhagic diathesis; Anticoagulation; Severe decompensation; Infections / inflammations (e.g. infections of the pleural cavity, florid skin infections in the area of the puncture site); Restricted access route to the pleural cavity (e.g. adhesions in the pleural space, malformation of the cardiac vessels, Altered anatomy of the chest wall); Emphysema of the lung

 **Intended user and patient target group:**

The product may only be inserted, applied and removed by Healthcare Professionals, i.e. by qualified physicians or by qualified medical personnel under the instruction of a qualified physician. Of particular importance are the procedures in anaesthesia, intensive care, internal medicine, pneumology (pulmonology), surgery and emergency medicine and for all other patients with the listed indications and where the user with the required qualifications considers this necessary. The medical device can be used on both adults and children without any fundamental distinction as to age, anatomy or physiology, taking into account the contraindications. The anatomical and physiological conditions of the patient must be checked by the attending Healthcare Professional before using the product.

 **General notes and precautions:**

- **The product is to be used sterile only!**
- If the product is delivered in a non-sterile state, the ***sterilization instructions*** must be taken into account (the sterilization instructions do not apply to products delivered in a sterile state in which the sterile state is no longer given (e.g. due to damage to the sterile barrier system or use of the product)).
- The product is intended for single use only and must not be resterilized. In case of reuse or reesterilization, it cannot be guaranteed that biological, physical, chemical and functional properties of the product are still given. Reuse or reesterilization may lead to malfunction of the product, which in turn may lead to injury, illness or death of the patient.
- The product may only be used under absolutely aseptic conditions.
- The product may only be used in accordance with the intended purpose described.
- The medical methods and procedures described in this manual do not reflect all accepted medical practices nor are they a substitute for the experience and judgment of the physician in treating his or her patients.
- Users must be aware of the risks, complications, and undesirable side effects associated with pleural puncture procedures, including but not limited to those listed in this instruction manual.
- In case of treating pregnant women, the use of the product must be checked by a doctor before application.
- Improper handling can damage the catheter (e.g. shearing of the catheter) with the consequence that parts of the catheter could remain in the patient and cause him/her harm.
- Do not use damaged products or products from damaged packaging, as this could result in the product no longer functioning properly, parts getting into the patient's body or the sterile barrier system no longer being intact with the consequence that sterility is no longer guaranteed.
- Do not use products that have exceeded the expiry date on the packaging (label).
- When using additional or other components that are not included in the packaging, care must be taken to ensure safe and firm connection, compatibility with the components contained in the packaging, as well as the instructions of the respective manufacturers.
- The connections may have become loose during transport and storage or during the sterilization process. Check all connections before use and tighten the connections if necessary.
- If a **malfunction** of the medical device occurs, the application must be stopped immediately and the manufacturer must be informed. Defective products must be kept and sent to the manufacturer for investigation.
- **Serious incidents** related to the medical device must be reported immediately to the manufacturer and the competent authorities, and the products must be kept and sent to the manufacturer for investigation.

 **Possible side effects:**

- Potential allergic reaction to materials with significant body contact: Catheter (polyethylene, barium sulphate); puncture needle (stainless steel).

 **Possible risks, complications and further notes:**

(a) The terms „adult“, „paediatric“ and „neonate“ in the specification description are for guidance. Which product (size) is most appropriate for the patient must be assessed by the clinician prior to use according to the patient's anatomical and physiological conditions; (b) Pleural drainage malposition: 1. Insertion too deep: The drain causes chest or shoulder pain when pressure is applied to the parietal pleura, and the drain is also compromised in its draining function; 2. Extrathoracic malposition: The drain has slipped past the rib and is now in subcutaneous or submuscular tissue. → This is noticeable by the fact that the fluid column in the drainage system no longer moves in breath synchronous; 3. Intrapulmonary malposition: This malposition is mainly possible in the case of adhesions and represents a lung injury; 4. Injury to intercostal

nerves and vessels which run along the lower edge of the respective rib. **Recommendation:** The thoracic drain should always be placed at the upper edge of the rib in order not to injure the nerves and vessels; (c) Before each drainage procedure, check the position of the catheter; (d) Sagging of the tubes / reduction of suction: Make sure that the tubes do not sag, as fluids or secretions collect in the loops (siphon), which reduces suction; (e) Infections / inflammations (e.g. at the puncture site, intrapleural wound infection); (f) skin emphysema; (g) injury to intra-abdominal or intrathoracic organs (visceral injuries) (e.g. heart, lungs, lung parenchyma, diaphragm), intrapleural injury or nerve damage; (h) haemorrhage; (i) circulatory distress when aspirating excessive amounts of pleural effusion; (j) pneumothorax if used incorrectly; (k) for pleural drainage, a suction of -10 to -20 cm H<sub>2</sub>O on the patient is sufficient for aspiration and to maintain intrapleural pressure. Unless clinically necessary, negative pressures greater than -50 cm H<sub>2</sub>O should not be applied to the patient. When used on children and neonates, a lower suction level may be necessary. It is recommended to apply a suction of -5 to -10 cm H<sub>2</sub>O to the patient in children and neonates. Excessive suction during drainage may result in severe pain, lung injury, injury to surrounding blood vessels, syncope or re-expansion oedema;

**(l) Alternative treatment options:**

- (1) For effusions < 100 ml in the pleural cavity, consider drug therapy as an alternative treatment option;
- (2) For effusions in the pleural cavity: Thick fluids may cause blockage of the catheter and puncture cannula. In this case, the use of drains with a larger diameter is recommended; (3) If there are appropriate indications, as well as recurrent accumulations of fluid or air in the pleural cavity in the case of prolonged drainage or after drainage has been completed, a thoracoscopy or thoracotomy may be useful to clarify the cause and for treatment. Alternatively, the use of a minithoracotomy in combination with upstream imaging for exploration of the pleural cavity can provide essential information about the intrathoracic findings, with the possibility of immediately placing a drain.

**Application steps:**

1. Open the packaging and remove the components (→ picture 1: possible package contents).
2. Carefully remove the protective tube from the puncture needle (→ picture 2).
3. Check the mobility between the catheter and the puncture needle by slightly turning the needle (→ picture 3).
4. Close the Three-Way-Stopcock in the catheter direction (off-position towards the metal tube connection from the Three-Way-Stopcock) (→ picture 4).
5. Determine the puncture site depending on the indication and the previously determined localisation of the medium to be aspirated (e.g. by imaging procedures).
6. Disinfect the puncture site and administer local anaesthesia.
-  **Attention: Use only alcohol or iodine-based disinfectants. Other disinfectants may weaken the material on contact with the catheter, which could affect the catheter properties. This may lead to leaks or air aspiration.**
7. Perform the puncture by rotational advancement. This eliminates the need for a stab incision (→ picture 5). Supportive ultrasound examination can be helpful for exact positioning.
8. On reaching the pleural space:
  - a) if the pleural pressure is negative, air may escape from the thrust sheath and the thrust sheath may visibly contract (→ picture 6).
  - b) if the pleural pressure is positive, the sheath may fill with air or effusion (→ picture 7).
9. Push the pleural drain provided with lateral eyes into the chest wall to the desired length, but at least until the double marking on the catheter in the puncture cannula disappears (→ picture 8 / 9).
10. Pull back the puncture needle together with the sterile sheath over the pleural drain fixed at the same time with the other hand. (→ picture 10).

 **Attention: Never pull the catheter back through the needle, as there is a risk of the catheter being sheared off by the needle or material of the catheter being abraded, thus causing the catheter to leak or plastic particles to enter the patient! (→ picture 11)**

11. Fix the catheter (e.g. with adhesive plaster) to the skin (→ picture 12).
-  **Attention: Immediately after placement of the catheter, it must be temporarily fixed so that it cannot be pulled into the pleural cavity by a possible negative pressure! The temporary fixation may only be released when the final fixation of the catheter is ensured (see point 15).**
12. Clamp the pleural drain with the thumb and forefinger of one hand and remove the distal plug with the other hand (→ picture 13).
13. Connect the Three-Way-Stopcock with the change-over tap pointing in the direction of the metal tube connection to be inserted (closed position of the catheter) (→ picture 14).
-  **Attention: The catheter must be pushed onto the metal mandrel as far as possible. Hold the catheter at the end to avoid kinking the tube. Afterwards, it is essential to check that the catheter is firmly seated on the metal mandrel (→ picture 15).**
14. A connection tube, a Three-Way-Stopcock, a syringe LL, a collection bag can be connected for the purpose of drainage, instillations or flushes. A Heimlich-valve PNEUMOVENT® can also be connected to prevent air from entering the thoracic cavity and thus further pneumothorax. In addition, there is also the possibility of connecting a manometer for pressure measurement.

 **Attention: Close the bottom drain valve of the collection bag before use! (→ picture 16)**

Different versions of the valve are possible.

15. Fix the catheter in a suitable place. The catheter must be fixed in such a way that it cannot detach from the connection on the patient's body (→ picture 17 / 18).

 **Attention: Only after this fixing the attachment from point 11 may be detached.**

16. If necessary, check the position of the catheter by performing imaging procedures.
17. Start the application.

Depending on the suction applied, the drain allows up to 0.2 litres of liquid per minute to be sucked out.
18. When the application is complete, release the fixation and hold the catheter firmly.

 **Attention: Be sure to hold the catheter firmly during and after detaching the fixation so that it cannot be pulled into the pleural space by any negative pressure! (→ picture 19)**

19. Carefully pull the catheter out of the patient and protect the puncture site immediately afterwards (e.g. with a gauze swab and bandage) (→ picture 20).

 **Attention: Immediate closure of the puncture site (e.g. with gauze swab and dressing) is essential to avoid a recurrence of a pneumothorax and infections.****Disposal:**


After use, the product must be disposed of in accordance with the regulations for infectious waste or in accordance with national or regional regulations. When disposing of sharp-edged products (e.g. puncture needle), there is a risk of cutting resulting in potential infections and microbiological hazards.

**Excluding of liability:**

Due to the biological differences of the patients, the effectiveness cannot be fully guaranteed. As we have no control over the conditions of use, diagnosis and indication as well as application and use of the product, we can neither guarantee success nor can we exclude the occurrence of side effects. **intra special catheters GmbH** accepts no liability for damage or injury of any kind resulting from improper use, reuse or modification of the product.

**Symbols:**

Not all of the following listed symbols are necessarily applicable to a product. The labeling of the product is authoritative. The symbols used on the packaging and label have the following meaning:

Further information on the product can be found at the contact details provided.

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