

CLEANTEX®

pouches

Sterilization: Compatible with steam sterilization process and ethylene oxide (EO) sterilization process.

Not suitable for Gamma or Plasma irradiation process.



DESCRIPTION :

Flat sterilization pouches.

One side 60gsm non woven - one side lilac tinted printed laminate made of 12µm polyester and 38µm polypropylene.

Preformed Sterile Barrier System according to EN ISO 11607-1.

Pouches codes: SCX***X***

- Natural rubber latex free
- TSE/BSE : compliant to EMEA/410/01

APPLICATION

Designed for the packing of bulky medical devices.
Used as primary packaging (SBS) and / or secondary packaging (protective packaging).

- To close the packaging and create a Sterile Barrier System, the use of a validated heat sealing machine is recommended according to EN ISO 11607-2.
- Make sure to keep a sufficient area to allow the passage of the sterilizing agent (fill to 2/3 max capacity within the packaging)
- Position the packaging plastic to plastic or porous to porous sides together in the sterilizer.
- Recommended sealing parameters: 170 – 185 °C
Notes: For special material/sealer equipment it might be necessary to vary from the limit values (contact the machine manufacturer and follow their IFU if necessary).
- At the opening:
 - Remove the attachment points
 - Peel chevron side

PERFORMANCES

- **Excellent mechanical strength**
- **Great flexibility and convenience in packaging**
- **Microbial barrier according to ISO 11607-1**

STORAGE

- It is recommended that the products are kept in the original, closed transport carton and are stored in dry and clean conditions protected from direct sunlight and excessive moisture.
- 5 years from the manufacturing date, provided the above storage conditions are met.

CONFORMITY

NF EN ISO 11607-1

ORIGIN

Made in France



CE MARK

Class I medical device (accessory) according to MDR 2017/745/EU

TECHNICAL DATA :

Film web : 12 µm lilac polyester laminated to a 38 µm "shatterless" polypropylene, printed intrafilm.

Porous web : 60 gsm non woven, tested and validated according to EN 868-2.

Process indicator for Steam sterilization according to ISO 11140-1 (class1).

CLEANTEX®

PROPERTIES	UNITS	METHODS	MIN	TYPIC	MAX
SUBSTANCE	gsm	ISO 536	57	60	63
HIGH POROSITY	l/min/1dm²	EN 868-2:1999	6	18	
THICKNESS (BULKY MATERIAL)	µm	ISO 12625-3		220	
ROUGHNESS FF	ml/min	ISO 8791-2	250	350	500
ROUGHNESS FT	ml/min	ISO 8791-2	250	350	500
TENSILE STRENGTH MD	kN/m	ISO 1924-2	1,60	2,30	
TENSILE STRENGTH CD	kN/m	ISO 1924-2	0,65	1,00	
WET TENSILE STRENGTH MD	kN/m	ISO 3781	1,2	2	
WET TENSILE STRENGTH CD	kN/m	ISO 3781	0.5	0,8	
BURST STRENGTH	kPa	ISO 2758	160	230	
WET BURST	kPa	ISO 3689	90	200	
TEARING STRENGTH MD	mN	ISO 1974	750	1150	
TEARING STRENGTH CD	mN	ISO 1974	1000	1700	
STRETCH MD	%	ISO 1924-2	5	7	
STRETCH CD	%	ISO 1924-2	7	10	
MASON JAR	min	WSP 80-11	75	130	
WATER REPELLENCY	s	EN 868-2:20017	20	30	
FLUORESCENCE	pts/dm²	EN 868-2:20017		0	

LAMINATE

PROPERTIES	UNITS	METHODS	VALUES
POLYESTER SUBSTANCE	gsm	SPS	16 to 17.4
POLYESTER THICKNESS	µm	SPS	12 ± 5 %
ADHESIVE SUBSTANCE	gsm	SPS	1.65 ± 0.15
POLYPROPYLENE SUBSTANCE	gsm	SPS	32.5 to 35.9
POLYPROPYLENE THICKNESS	µm	SPS	36 to 40

POUCH

PROPERTIES	UNITS	METHODS	VALUES
SEAL STRENGTH	N/15mm	EN 868-5	≥3
DIMENSIONS	Mm	SPS	Nominal dimensions ±3