

Oxygen Masks with oxygen reservoir



REF: 01.000.08.131, 01.000.08.130

Product group: O2 Therapy Accessories

Intended use:

The oxygen masks are used to administer high concentrations of oxygen without rebreathing.

Product description:

The oxygen masks comprise a mask body with elastic strap, a reservoir bag (600 ml), a clip that adjusts to the bridge of the nose, and an oxygen supply tube with standard connectors and 2 valves.

Sterile: No

REF	Application area	Length (oxygen tube)	Selling unit quantity	Selling unit GTIN	Packing unit quantity	Packing unit GTIN
01.000.08.131	Child	210 cm	50 pcs	04034342119959	50 pcs	04034342119959
01.000.08.130	Adult	210 cm	50 pcs	04034342010508	50 pcs	04034342010508

Materials

Product	PVC, PP, aluminium, silicone
Packaging	PE bag, corrugated cardboard
Latex	No
DEHP/BBP/DBP	No

Additional information

Regulatory basis	Directive 93/42/EEC
Notified Body	TÜV SÜD Product Service GmbH (0123)
Classification	Ila
UMDNS	12-449
eClass	34241203
Manufacturer	P. J. Dahlhausen & Co. GmbH
Single-use Product	Yes
Sterilisation	-
Shelf life	5 Years
Instructions for Use	Yes