

**biliary drainage bag without secretion tube**

REF: 05.000.22.704

Unit: 50 pcs.

Packing Unit: 50 pcs.

**Description**

The biliary drainage bag with drainage connector, drain valve and fixation hook, is for direct connection to a drainage system in sizes Ch 12, 14, 15 and 16. The bag has a filling capacity of 500 ml.

**Material Composition**

<u>Product:</u>		PVC, PE
<u>Packing:</u>	Individual Packing:	PE Bag
	Outer case	Corrugated cardboard

**Sterilisation**

EO-Gas