

Umbilical Cord Clamps



REF: 07.044.00.000

Product group: Sterile Umbilical Cord Clamps

Intended use:

The umbilical cord clamp is used to securely close the postnatally severed umbilical cord in newborns. The umbilical cord clamp remains in place until the umbilical cord detaches by itself.

Product description:

The umbilical cord clamp consists of two toothed branches. The branches are joined together on one side. The open side features a lock which cannot be opened without aids once it has been audibly clicked into place.

Sterile: Yes

REF	Selling unit quantity	Selling unit GTIN	Packing unit quantity	Packing unit GTIN
07.044.00.000	100 pcs	04034342046316	600 pcs	04034342367633

Materials

Product	POM
Packaging	PP/PE, Paper, PE, corrugated cardboard
Latex	No
DEHP/BBP/DBP	No

Additional information

Regulatory basis	Directive 93/42/EEC
Notified Body	TÜV SÜD Product Service GmbH (0123)
Classification	Is
UMDNS	10-876
eClass	34280104
Manufacturer	P. J. Dahlhausen & Co. GmbH
Single-use Product	Yes
Sterilisation	EO gas
Shelf life	3 Years
Instructions for Use	Yes