

Chest Aspiration Set

REF: 07.048.00.000

Unit: 10 pcs.

Packing Unit: 30 pcs.

Description

The Chest Aspiration Set enables the puncture of Pleura for secretion drainage. The set is consisting of: secretion bag 2 l, tube with 3 way connector, puncture needle and syringe.

Material Composition

<u>Product:</u>	tube:	PVC
	syringe:	PP
	puncture needle:	stainless steel, ABS, PP, PE size: 2,1 x 80 mm (14G x 3 1/8)
	secretion bag 2l:	PVC
	3 way connector:	ABS, PC
<u>Packing:</u>	Individual Packing:	Paper, foil polyamide / polyethylene
	Outer case:	Corrugated cardboard

Sterilisation

EO-Gas