



Product Main Number 1000031138

2 layers, 210x80cm, white, airlaid/LDPE



### Specifications

Color	White
Single or multiple use	Single use
Material	Airlaid, LDPE
Length/depth	210 cm
Width	80 cm
Layer quantity	2 layers
Certificates	CE, MD.
CE Class (Medical Devices)	Class I
Directives, regulations and acts	MDR (EU) 2017/745
Storage Instructions	Store dry, room temperature and no direct sunlight.
Product Disposal Instructions	Can be disposed of with normal household waste sorted according to local regulations.
Packaging Disposal Instructions	Can be recycled or incinerated.

### Packaging data

Unit	Contains	Length	Width	Height	EAN
c/l	3 pck	54 cm	28 cm	24 cm	7331246006604
pck	15 pcs				7331246006598



Regulation (EU) 2017/745 is a regulation of the European Union on the clinical investigation and sale of medical devices for human use. It repeals Directive 93/42/EEC, which concerns medical devices, and Directive 90/385/EEC, which concerns active implantable medical devices, on 26 May 2021.



The CE mark guarantees that a product is safe to use and complies with all safety precautions. CE stands for Conformité Européenne (European Conformity) and is mainly found on electronic equipment, safety equipment, construction products, and medical devices.



Your guarantee that our products have been tested for ingredients which can potentially cause skin irritation. Tests conducted by an independent third-party laboratory.



The product is only for single-use.



Barriers lengthwise on the product can help reduce incontinence related leakages onto beds and other furniture or surfaces.