



17601

Short touch screen cover

SIZE: **53 x 45 cm**

QUANTITY: **80**

CONTAINS LATEX: **No**

STERILE: **Yes**

STERILIZATION METHOD: **EO**

MATERIAL: **Polyethylene**

GTIN PRIMARY: **50748426084901**

GTIN SECONDARY: **00748426084906**

GTIN SALES: **50748426084901**

COMMODITY CODE: **3926909790**

INTENDED USE: **covering equipment**

SHELF LIFE: **5 Years**

DEHP/PHTHALATES: **No**

GMDN CODE: **12535**

CE NOTIFYING BODY: **CE0482**

MEDICAL CLASSIFICATION: **Is**

BOX LENGTH (MM): **425**

BOX WIDTH (MM): **320**

BOX HEIGHT: **342**

BOX M3: **0.046512**

BOX WEIGHT (GRAM): **17300**

LABEL LANGUAGES: **BG,CS,DA,DE,EL,EN,ES,ET,FI,FR,HR,HU,IS,IT,LT,LV,NL,NO,PL,PT,RO,RU,SK,SL,SR,SV,TR**

2X STERILISATION ALLOWED?: **No**

ADDRESS OF LEGAL MANUFACTURER: **Sorbonne Centre F20, Mosta, Technopark, Mosta, MST3000 Malta**

BISPHENOL A: **No**

IFU: **Yes**

SINGLE USE: **Yes**

CASES PER PALLET CONVERSION: **18**

INFORMATION LAST UPDATED ON: **21 Jun 2022**
