

ECG Electrode

REF: 19.000.00.400

Unit: 30 pcs

Packing Unit: 1.500 pcs.

Description

Foam ECG electrode x-ray negative, oval, 40mm diameter.

The electrode is MRT suitable (MRI safe up to 3 Tesla)

The electrodes are free of latex and tested according to ISO-10993, biocompatible

Material Composition

<u>Product:</u>	carrier:	PE-foam with med. Acryl adhesive (40 g/m ²)
	sensor rivet:	carbon, Ag/AgCl coating
	press-stud:	carbon
	gel:	Hydrogel for ECG
	cover:	siliconized PET

<u>Packing:</u>	pouch:	paper/Alu/PE (steamtight)
	box:	corrugated cardboard
	outer case:	corrugated cardboard

Sterilisation

Non-sterile