



**Alleset**  
4142 Industry Way  
Flowery Branch, GA 30542  
**Main:** (678) 866-0550  
**Fax:** (678) 866-0551

**To Whom It May Concern**

**From:** Alleset Inc.

**Date:** May 10th, 2024

**RE:** latex free

We (GRI) hereby declare that GRI produced product listed below can be placed in a Warming device set at 170 Deg F (76 deg C) for a period of 8 hours without any noticeable change in its appearance.

It is not recommended to heat the blanket with temperatures greater than 170 deg F or to put them through repeated heating process.

| <b>PRODUCT CODE</b> | <b>Product Description</b>           | <b>STERILITY</b> |
|---------------------|--------------------------------------|------------------|
| <b>68-4115175</b>   | Standard Warming Blanket_130cmX190cm | Non-Sterile      |

If you have any questions, please don't hesitate to contact Wendy Liu at [wliu@gri-allese.com](mailto:wliu@gri-allese.com).

Sincerely,

A handwritten signature in black ink that reads "Wendy Liu". The signature is written in a cursive, flowing style.

Wendy Liu  
Global Product Manager  
E: [wliu@gri-allese.com](mailto:wliu@gri-allese.com)