

TECHNICAL DATA SHEET

Description

Reference: 69601 XXL

Impervious Film Laminate- Blue, Disposable Fluid-Resistant, Potential for Soaking and Heavy Fluid Contact in Front, Perforated Open Back.

Dimensions

L XXL

Length (from shoulder to hem) 121 cm 121 cm

Sleeve length (underarm to cuff/hook) 61 cm 62 cm

Neck length 27 cm 27 cm

Body width (chest) 71 cm 84 cm

Width (overall) : 119 cm 140 cm

Properties

Basis weight of base material: 38 g/m²

Material was tested as per ASTM F 1671- Viral Penetration and as per ASTM F 739-99a- Resistance to Permeation by Liquids and Gases.



Indication

To prevent the transfer of infective agents between the operating room personnel and the patient.
Single use product.

Counter Indication

none in particular.

Main Materials

Base material: Spunbond-Film Laminate (polypropylene)

Cuffs and Threads: Polyester

Neck Closure: Fastening Tape with Release Liner.

Does not contain Natural rubber latex.

Does not contain DEHP.

Sterilisation

non-sterile

Shelf Life

none specified

Packaging

69600, 69601, 69602, and 69603: Shipping case of 100 units.
10 units are placed within 1 polyethylene bag,
10 polyethylene bags are placed within 1 shipping case.
Bar coding: GS1-128 symbology, linear, on shipping case and bags.

Manufacturing

Products are manufactured in China.

The quality system of the manufacturing sites is ISO 13485 compliant.

Regulatory Information

Product CE marked as per Regulation (EU) 2017/745 on Medical Devices.

Class of the device: I.

Storage

Store in a dry and cool place, away from intense sources of heat.
Keep as much as practicably possible in its shipper box.