

TECHNICAL DATA SHEET

Description

Reference: 69766

3-layers SMS Fabric, Sleeveless, Full Length with Overlapping Back Panel Secured by Waist Belt, Blue

Dimensions: Length (front, neck to bottom) 117cm Neck length 18 cm Body width (chest) 66cm Opening for arm 29cm Belt length 168cm

Properties: Basis weight of base material: 28g/m²

Indication

Gown ideal for patient comfort and privacy. Single use product.

Counter Indication

None in particular.

Main Materials

Base material: SMS Fabric (polypropylene), topically treated for alcohol repellency and static dissipation. Does not contain Natural rubber latex. Does not contain DEHP.

Sterilisation

non-sterile

Shelf Life

None specified.

Packaging

Shipping case of 100 units. 10 units are placed within 1 polyethylene bag and 10 polyethylene bags are placed within the shipping case. Bar coding: GS1-128 symbology, linear on shipping case and bags.

Manufacturing

Product is manufactured in Mexico or USA. The quality system of the manufacturing site is ISO13485 compliant.

Regulatory Information

Product CE marked as per Regulation (EU) 2017/745 on Medical Devices.

Class of the device: I.

Storage

Store in a dry and cool place, away from intense sources of heat. Keep as much as practicably possible in its shipper box.