



8888MEU

C-arm drape for Fluoroscans

SIZE: **137 x 198 cm**

QUANTITY: **40**

CONTAINS LATEX: **No**

STERILE: **Yes**

STERILIZATION METHOD: **EO**

MATERIAL: **Polyethylene**

GTIN PRIMARY: **50748426099196**

GTIN SECONDARY: **00748426099191**

GTIN SALES: **50748426099196**

COMMODITY CODE: **3926909790**

SHELF LIFE: **5 Years**

DEHP/PHTHALATES: **No**

GMDN CODE: **12535**

CE NOTIFYING BODY: **CE0044**

MEDICAL CLASSIFICATION: **Is**

BOX LENGTH (MM): **540**

BOX WIDTH (MM): **440**

BOX HEIGHT: **285**

BOX M3: **0.067716**

BOX WEIGHT (GRAM): **11350**

COUNTRY OF ORIGIN: **Dominican Republic**

LABEL LANGUAGES: **CS,DA,DE,EL,EN,ES,FI,FR,HU,IT,NL,PL,PT,SV**

2X STERILISATION ALLOWED?: **Yes**

ADDRESS OF LEGAL MANUFACTURER: **1 Ecolab Place, St. Paul, MN 55102 U.S.A.**

BISPHENOL A: **No**

IFU: **Yes**

PVC: **No**

REACH COMPLIANCE: **Yes**

SINGLE USE: **Yes**

CASES PER PALLET CONVERSION: **14**

INFORMATION LAST UPDATED ON: **21 Jun 2022**
