



8911

Towel drape / Keyboard cover

SIZE: **51 x 61 cm**

QUANTITY: **40**

CONTAINS LATEX: **No**

STERILE: **Yes**

STERILIZATION METHOD: **GAMMA**

MATERIAL: **Polyethylene**

FOLDED/ROLLED: **Folded**

GTIN PRIMARY: **50748426023511**

GTIN SECONDARY: **00748426023516**

GTIN SALES: **50748426023511**

COMMODITY CODE: **3926909790**

SHELF LIFE: **4 Years**

GMDN CODE: **47783**

CE NOTIFYING BODY: **NA**

MEDICAL CLASSIFICATION: **NA**

BOX LENGTH (MM): **295**

BOX WIDTH (MM): **230**

BOX HEIGHT: **145**

BOX M3: **0.00983825**

BOX WEIGHT (GRAM): **1200**

LABEL LANGUAGES: **CS,DA,DE,EL,EN,ES,FI,FR,HU,IT,NL,PL,PT,SV**

2X STERILISATION ALLOWED?: **No**

ADDRESS OF LEGAL MANUFACTURER: **1 Ecolab Place, St. Paul, MN 55102 U.S.A.**

BISPHENOL A: **No**

IFU: **Yes**

REACH COMPLIANCE: **Yes**

SINGLE USE: **Yes**

CASES PER PALLET CONVERSION: **48**

INFORMATION LAST UPDATED ON: **21 Jun 2022**
