

Secretion Bag 1,5 l

REF: 94.000.00.001

Unit: 1 pc

Packing Unit: 100 pcs

Description

Sterile secretion bag (universal) with connector. The bag enables the sampling of body fluids out of a catheter and/or drainage.

Length of the tube: 86 cm

Material Composition

Product: bag: PVC (DEHP)
connector: ABS, PS, PE, Polyisoprene

Packing: individual packing: PE-bag, paper, foil PE/PA
outer case: corrugated cardboard

Sterilisation

EO-Gas