



APC1290

Adhesion General Purpose Probe Cover with Gel

SIZE: **15 x 122 cm**

WIDTH IN CM: **15**

LENGTH IN CM: **122**

QUANTITY: **20**

CONTAINS LATEX: **No**

STERILE: **Yes**

STERILIZATION METHOD: **ETO**

MATERIAL: **Ethyl Methyl Acrylate/anti-stat blend**

FOLDED/ROLLED: **Folded**

GTIN PRIMARY: **50748426107372**

GTIN SECONDARY: **30748426107378**

GTIN SALES: **50748426107372**

COMMODITY CODE: **3926909790**

SHELF LIFE: **The expiration date must match the expiration date on the included Gel Sachet. Consult QA for detailed information.**

DEHP/PHTALATES: **No**

GMDN CODE: **44713**

CE NOTIFYING BODY: **CE0044**

MEDICAL CLASSIFICATION: **IIa**

BOX LENGTH (MM): **286**

BOX WIDTH (MM): 222

BOX HEIGHT: 140

BOX M3: 0,008888888

BOX WEIGHT (GRAM): 1500

LEGAL MANUFACTURER: Microtek Medical, Inc.

COUNTRY OF ORIGIN: Dominican Republic

LABEL LANGUAGES: CZ, DE, DK, EN, ES, FR, GR, HU, IT, NL, PL, PT, RU, SE

2X STERILISATION ALLOWED?: Yes

ADDRESS OF LEGAL MANUFACTURER: 602 Lehmberg Road, Columbus, MS 39702 U.S.A.

BISPHENOL A: No

CE CERTIFIED: Yes

CE CERTIFICATE: CE0044

IFU: Yes

PVC: No

REACH COMPLIANCE: Yes

SINGLE USE: Yes

SINGLE PACKED: Yes