

Oxygen Mask

Medium Concentration with standard connector

Intended use:

Medium concentration masks are intended for administration of oxygen when precise control of inspired oxygen concentration is not mandatory. During inspiration the patient draws gases both from oxygen flowing into the mask through the tubing, as well as from room air via ports on the sides of the mask. Oxygen concentration delivered is 40%-60% depending on the patient's breathing. The masks are connected directly to a compressed air or oxygen supply. The products are intended for single use only.



Product range:

	REF*	Item	Size	Description	Tube length
 9-102MM  9-222MM  9-262MM	9-102MM 9-102MMV (Vanilla scented mask)	Medium Concentration Oxygen Mask with tubing	Adult	Medium Concentration Oxygen Mask, Elongated mask with nose clip and elastic strap, 2.1m (7') No-crush oxygen tube, Standard connector	2.1m (7')
	9-222MM 9-222MMV (Vanilla scented mask)	Medium Concentration Oxygen Mask with tubing	Pediatric	Medium Concentration Oxygen Mask, Elongated mask with nose clip and elastic strap, 2.1m (7') No-crush oxygen tube, Standard connector	2.1m (7')
	9-262MM	Medium Concentration Oxygen Mask with tubing	Infant	Medium Concentration Oxygen Mask, Mask with nose clip and elastic strap, 2.1m (7') No-crush oxygen tube, Standard connector	2.1m (7')
	9-1023MMV (Vanilla scented mask)	Medium Concentration Oxygen Mask with tubing	Adult	Medium Concentration Oxygen Mask, Elongated mask with nose clip and elastic strap, 3m (10') No-crush oxygen tube, Standard connector	3m (10')

* Products may not be available in all countries. Please contact your regional ConvaTec sales representative for further information.

Material specifications:

Contains latex? No
Contains phthalates? No
Contains bisphenols? No
Mask Polyvinyl chloride
Elastic strap Polyester
Nose clip Aluminum
Oxygen nipple connector High density polyethylene
Oxygen tube Polyvinyl chloride
Tube connector Polyvinyl chloride

Certifications:

CE-marking? Yes (CE 0120)
Product classification Class IIA (EU)
Quality system certified to EN ISO 13485:2012

Sterilization:

Sterile? No
Shelf life 5 years

Oxygen Mask

Medium Concentration with standard connector



Packaging specifications:

REF	Primary package	Secondary package/Polybag	Tertiary package/Shipper			Pallet
	Qty, pcs	Qty, pcs	Qty, pcs	Dimensions, mm (LxWxD)	Gross Wt., kg	Qty, pcs
9-102MM	1	10	50	400 x 380 x 260	4.00	900
9-102MMV	1	10	50	400 x 380 x 260	4.00	900
9-222MM	1	10	50	400 x 380 x 260	3.28	900
9-222MMV	1	10	50	400 x 380 x 260	3.28	900
9-262MM	1	10	50	400 x 300 x 240	3.05	1,600
9-1023MMV	1	10	50	400 x 380 x 260	4.55	900

EAN codes:

REF	Primary package	Secondary package	Tertiary package
9-102MM	(01)06945099301635	(01)26945099301639	(01)46945099301633
9-102MMV	(01)06945099301857	(01)26945099301851	(01)46945099301855
9-222MM	(01)06945099301642	(01)26945099301646	(01)46945099301640
9-222MMV	(01)06945099301864	(01)26945099301868	(01)46945099301862
9-262MM	(01)06945099301871	(01)26945099301875	(01)46945099301879
9-1023MMV	(01)06945099301888	(01)26945099301882	(01)46945099301886

Packaging material:

Primary package	Poly bag - Low density polyethylene / Insert - 80g paper
Secondary package	Poly bag - Low density polyethylene / Sticker - Coated paper
Tertiary package	Carton - Corrugated board / Sticker - Coated paper

Storage and handling:

	FRAGILE - Contents of the transport package are fragile and shall therefore be handled with care
	KEEP DRY
	THIS WAY UP - Indicates correct upright position of the transport package
	DO NOT USE HAND HOOKS - Hooks are prohibited for handling the transport package



Ningbo Shengyurui Medical Appliances Co., Ltd.

No. 138 Binhaisi Rd.
Hangzhou Bay New Zone, 315336, Ningbo City
Zhejiang Province, P.R. China

Made in China



Shanghai Intl. Holding Corp. GmbH

Eiffestrasse 80, 20537 Hamburg
Germany

Distributed by



ConvaTec Limited
First Avenue, Deeside Industrial Park
Deeside, Flintshire, CH5 2NU, UK

www.convatec.com

For further information on these products, please contact your local ConvaTec sales representative.

With reservations for misprints and alterations in product specifications. Printed version not controlled.

All trade marks and logos are the property of a ConvaTec group company.

© 2018 ConvaTec Inc. AP-018908-MM, ver. JUL 2018