

Disposable Prep Razor

REF:	10.000.00.010	double edged
	10.000.00.011	single edged
Unit:	200 pcs.	
Packing Unit:	200 pcs.	

Description

Disposable Prep Razor, sterile, double or single edged. The razors enable the removal of body hair.

Material Composition

<u>Product:</u>	handle:	PS
	blade:	stainless steel
	protection cap:	PP

<u>Packing:</u>	individual packing:	paper, foil PA/PE
	box:	cardboard
	case:	corrugated cardboard

Sterilisation

EO-Gas