

## Technical Datasheet - IrrigaCare enema bag, 1.750l, 150 cm tube, non-sterile

### Description & Purpose:

- The IrrigaCare enema bag is a single-use product for the rectal administration of water or non-medicated enema solutions.
- Administration is performed via a rectal tube inserted into the colon.
- The product is intended to support the evacuation of the distal colon.
- It is intended for patients requiring an enema, for example in cases of constipation.
- The maximum continuous application time is 60 minutes.



### Overview of the product range:

Art. Number	Volume	Tube length	PZN
07/01/1750/150/GHC	1750 ml	150 cm	20367719

### Material / product details

Contains Latex	No
Contains DEHP	Yes
For single use only	Yes
Medical product	Yes

### Certificates

CE marking	CE
Product classification	I

### Sterilization

Sterile	No
Re-sterilizable	No
Durability	5 Years

### Packaging

Art. Nr.	Quantity / PU	Weight / PU	Packed	Purchase quantity
07/01/1750/150/GHC	200 Stück	8,5 kg	einzel steril verpackt	1 VE

### Storage and handling

Avoid excessive heat
Protect from frost
Protect from moisture
Single-use product

**GHC** German Health Care GmbH

Manufacturer: Conod Medical Co., Limited | No. 11  
Hongfeng Road, Baimao Industrial Park, Guli Town |  
215532 Changshu City, Jiangsu Province | PEOPLE'S  
REPUBLIC OF CHINA

EU-Rep.: Shanghai International Holding Corp. GmbH  
Eiffestr. 80, 20537 Hamburg - Germany  
CH-Rep.: The Tao of Excellence GmbH | Vorstadt  
26 | CH-8200 Schaffhausen - Switzerland

Importer: GHC German Health Care GmbH |  
Birkstück 15 | 24999 Wees | Germany |  
Tel. +49 (0) 4631 57 200-0 |  
info@ghc-medical.de | www.ghc-medical.de

Date 25.03.2026

Revision 001/03.2026