


# Intermittent Catheters



### Area of application

Intermittent Catheters are used to provide an intermittent pathway for short term drainage of urine from the bladder of male (Nelaton and Tiemann) and female (Female) patients. The products are for single use only.

### Product range

Size	Nelaton Catheters		Nelaton Catheters with Sheath		Tiemann Catheters	Female catheters
	Standard (Shore A 76)	SOFT (Shore A 64)	Standard (Shore A 76)	SOFT (Shore A 64)	Standard (Shore A 76)	Standard (Shore A 76)
CH/Fr.	Effective Length 370mm*	Effective Length 370mm*	Effective Length 370mm*	Effective Length 370mm*	Effective Length 370mm*	Effective Length 200mm*
6	11300118		11400118		11700118	11200118
8	11300218		11400218		11700218	11200218
10	11300318	11600318	11400318	11500318	11700318	11200318
12	11300418	11600418	11400418	11500418	11700418	11200418
14	11300518	11600518	11400518	11500518	11700518	11200518
16	11300618	11600618	11400618	11500618	11700618	11200618
18	11300718	11600718	11400718	11500718	11700718	11200718
20	11300818		11400818		11700818	

\*Length from tip to connector, Length of connector: 40mm

### Sterilization

Sterile, Sterilized by ethylene oxide, Shelf life: 5 years

### Certifications

<b>CE-marking</b>	CE 0123
<b>Product classification</b>	Class IIA (EU)
<b>Quality system certified to</b>	ISO 13485 Medical Device Single Audit Program (MDSAP)
<b>European Medical Device Nomenclature (EMDN)</b>	U01010502 - NELATON CATHETERS, NON SELF-LUBRICATING
<b>Global Medical Device Nomenclature (GMDN)</b>	34930

### Instructions for use

Included?	Yes
-----------	-----

### Material specifications

<b>Contains latex?</b>	No
<b>Contains phthalates?</b>	No
<b>Contains cotton?</b>	No
<b>Contains bisphenols?</b>	No
<b>Contains paper?</b>	Packaging only
<b>Catheter shaft and funnel</b>	Polyvinyl chloride
<b>Sheath</b>	Polyethylene



# Intermittent Catheters

Packaging information									
Reference	Primary pack./ Single peelpack		Secondary pack./ Shelf box		Tertiary pack./ Shipper			Gross weight kg	Pallet Qty. pcs
	Qty. pcs	Dimensions mm (LxW)	Qty. pcs	Dimensions mm (LxWxD)	Qty. pcs	Dimensions mm (LxWxD)	Volume m <sup>3</sup>		
11300118	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	5.40	12,000
11300218	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	5.60	12,000
11300318	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	6.40	12,000
11300418	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	6.92	12,000
11300518	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	7.70	12,000
11300618	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	7.15	12,000
11300718	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	8.17	12,000
11300818	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	10.10	12,000
11600318	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.50	12,000
11600418	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	6.92	12,000
11600518	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	7.70	12,000
11600618	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	7.15	12,000
11600718	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	8.17	12,000
11400118	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.50	12,000
11400218	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.50	12,000
11400318	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.50	12,000
11400418	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	6.92	12,000
11400518	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	7.70	12,000
11400618	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	7.15	12,000
11400718	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.50	12,000
11400818	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.50	12,000
11500318	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.50	12,000
11500418	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	6.92	12,000
11500518	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	7.70	12,000
11500618	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	7.15	12,000
11500718	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	9.50	12,000
11700118	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	5.60	12,000
11700218	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	6.00	12,000
11700318	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	6.40	12,000
11700418	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	7.00	12,000
11700518	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	7.60	12,000
11700618	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	7.20	12,000
11700718	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	8.23	12,000
11700818	1	485 × 33.5	100	500 × 175 × 120	600	520 × 375 × 370	0.0722	10.00	12,000
11200118	1	322 × 33.5	100	325 × 175 × 120	600	380 × 345 × 370	0.0485	4.20	18,000
11200218	1	322 × 33.5	100	325 × 175 × 120	600	380 × 345 × 370	0.0485	4.50	18,000
11200318	1	322 × 33.5	100	325 × 175 × 120	600	380 × 345 × 370	0.0485	4.70	18,000
11200418	1	322 × 33.5	100	325 × 175 × 120	600	380 × 345 × 370	0.0485	4.16	18,000
11200518	1	322 × 33.5	100	325 × 175 × 120	600	380 × 345 × 370	0.0485	5.30	18,000
11200618	1	322 × 33.5	100	325 × 175 × 120	600	380 × 345 × 370	0.0485	5.60	18,000
11200718	1	322 × 33.5	100	325 × 175 × 120	600	380 × 345 × 370	0.0485	6.20	18,000



# Intermittent Catheters

GTIN codes (for production until 30 June 2025)			
REF	Primary package	Secondary package	Tertiary package
11300118	06936881134517	16936881134514	26936881134511
11300218	06936881134524	16936881134521	26936881134528
11300318	06936881134531	16936881134538	26936881134535
11300418	06936881134548	16936881134545	26936881134542
11300518	06936881134555	16936881134552	26936881134559
11300618	06936881134562	16936881134569	26936881134566
11300718	06936881134579	16936881134576	26936881134573
11300818	06936881134586	16936881134583	26936881134580
11600318	06936881152979	16936881152976	26936881152973
11600418	06936881152986	16936881152983	26936881152980
11600518	06936881152993	16936881152990	26936881152997
11600618	06936881153006	16936881153003	26936881153000
11600718	06936881153013	16936881153010	26936881153017
11400118	06936881152849	16936881152846	26936881152843
11400218	06936881152856	16936881152853	26936881152850
11400318	06936881152863	16936881152860	26936881152867
11400418	06936881152870	16936881152877	26936881152874
11400518	06936881152887	16936881152884	26936881152881
11400618	06936881152894	16936881152891	26936881152898
11400718	06936881152900	16936881152907	26936881152904
11400818	06936881152917	16936881152914	26936881152911
11500318	06936881152924	16936881152921	26936881152928
11500418	06936881152931	16936881152938	26936881152935
11500518	06936881152948	16936881152945	26936881152942
11500618	06936881152955	16936881152952	26936881152959
11500718	06936881152962	16936881152969	26936881152966
11700118	06936881134616	16936881134613	26936881134610
11700218	06936881134623	16936881134620	26936881134627
11700318	06936881134630	16936881134637	26936881134634
11700418	06936881134647	16936881134644	26936881134641
11700518	06936881134654	16936881134651	26936881134658
11700618	06936881134661	16936881134668	26936881134665
11700718	06936881134678	16936881134675	26936881134672
11700818	06936881134685	16936881134682	26936881134689
11200118	06936881134418	16936881134415	26936881134412
11200218	06936881134425	16936881134422	26936881134429
11200318	06936881134432	16936881134439	26936881134436
11200418	06936881134449	16936881134446	26936881134443
11200518	06936881134456	16936881134453	26936881134450
11200618	06936881134463	16936881134460	26936881134467
11200718	06936881134470	16936881134477	26936881134474

# Intermittent Catheters



GTIN codes (for production from 1 July 2025 onwards)			
REF	Primary package	Secondary package	Tertiary package
11300118	06925816648892	16925816648899	26925816648896
11300218	06925816648908	16925816648905	26925816648902
11300318	06925816648915	16925816648912	26925816648919
11300418	06925816648922	16925816648929	26925816648926
11300518	06925816648939	16925816648936	26925816648933
11300618	06925816648946	16925816648943	26925816648940
11300718	06925816648953	16925816648950	26925816648957
11300818	06925816648960	16925816648967	26925816648964
11600318	06925816649103	16925816649100	26925816649107
11600418	06925816649110	16925816649117	26925816649114
11600518	06925816649127	16925816649124	26925816649121
11600618	06925816649134	16925816649131	26925816649138
11600718	06925816649141	16925816649148	26925816649145
11400118	06925816648977	16925816648974	26925816648971
11400218	06925816648984	16925816648981	26925816648988
11400318	06925816648991	16925816648998	26925816648995
11400418	06925816649004	16925816649001	26925816649008
11400518	06925816649011	16925816649018	26925816649015
11400618	06925816649028	16925816649025	26925816649022
11400718	06925816649035	16925816649032	26925816649039
11400818	06925816649042	16925816649049	26925816649046
11500318	06925816649059	16925816649056	26925816649053
11500418	06925816649066	16925816649063	26925816649060
11500518	06925816649073	16925816649070	26925816649077
11500618	06925816649080	16925816649087	26925816649084
11500718	06925816649097	16925816649094	26925816649091
11700118	06925816649158	16925816649155	26925816649152
11700218	06925816649165	16925816649162	26925816649169
11700318	06925816649172	16925816649179	26925816649176
11700418	06925816649189	16925816649186	26925816649183
11700518	06925816649196	16925816649193	26925816649190
11700618	06925816649202	16925816649209	26925816649206
11700718	06925816649219	16925816649216	26925816649213
11700818	06925816649226	16925816649223	26925816649220
11200118	06925816648823	16925816648820	26925816648827
11200218	06925816648830	16925816648837	26925816648834
11200318	06925816648847	16925816648844	26925816648841
11200418	06925816648854	16925816648851	26925816648858
11200518	06925816648861	16925816648868	26925816648865
11200618	06925816648878	16925816648875	26925816648872
11200718	06925816648885	16925816648882	26925816648889

# Intermittent Catheters



Packaging material	
Primary package	Coated paper, Film
Secondary package	Ivory Cardboard
Tertiary package	Double-wall corrugated fiberboard

Storage and handling	
	This way up
	Keep dry
	Fragile, handle with care
	Keep away from sunlight
	Do not re-sterilize
	Do not re-use
	Do not use if package is damaged and consult instructions for use
	Temperature limit (-15° to 49°C)
	Consult instructions for use or consult electronic instructions for use

**DISTRIBUTED BY**  
**UNOQUIP GMBH**  
 Schützengraben 20  
 8200 Schaffhausen  
 Switzerland  
[info@unoquip.com](mailto:info@unoquip.com)  
[www.unoquip.com](http://www.unoquip.com)

**SPONSORED IN AUSTRALIA BY**  
 UNOQUIP Pty Ltd  
 PO Box 1117  
 Windsor VIC 3181  
 1300 55 80 79  
[info@unoquip.au](mailto:info@unoquip.au)



**Certified Quality Management System**  
 DIN EN ISO 13485:2021  
 EN ISO 13485:2016 / A11:2021  
 ISO 13485:2016

Well Lead Medical Co., Ltd.  
 C-4 Jinhua Industrial Estate  
 Hualong, Panyu Guangzhou  
 511434 P.R. China

**EU REP** Shanghai Intl. Holding  
 Corp. GmbH  
 Eiffestrasse 80  
 20537 Hamburg  
 Germany

**CH REP** QS Engineering AG  
 Rütihof 1  
 CH 8820 Wädenswil  
 Switzerland

**UK Responsible Person (UKRP)** SUNGO Certification  
 Company Limited  
 3rd Floor, 70 Gracechurch  
 Street, London, EC3V 274HR  
 England, United Kingdom