




Home Care Leg Bags

Sterile Urine Drainage Bags



Area of application

Urine drainage bags are used for the collection and monitoring of urine output. The scale is approximate and for indication purposes only. The product is intended for single use only.

Product range								
Sterile Home Care Leg Bag with Bottom Outlet and Soft Backing								
Capacity, ml	Tube Length cm	Cut to fit	Needle-free Sample Port	Bottom Outlet*	Leg Staps*	Soft Backing	Packaging	Reference
350	30	-	Yes	Cross Valve	Yes	Yes	1/10/80	23510118
350	9	-	Yes	Cross Valve	Yes	Yes	1/10/80	23510218
500	9	-	Yes	Cross Valve	Yes	-	1/10/80	25010118
500	30	-	Yes	Cross Valve	Yes	-	1/10/80	25010218
500	30	-	Yes	Cross Valve	Yes	Yes	1/10/80	25010318
500	9	-	Yes	Cross Valve	Yes	Yes	1/10/80	25010418
500	30	-	Yes	Lever tap	Yes	-	1/10/80	25010518
500	9	-	Yes	Lever tap	Yes	-	1/10/80	25010618
500	30	-	Yes	Lever tap	Yes	Yes	1/10/80	25010718
500	9	-	Yes	Lever tap	Yes	Yes	1/10/80	25010818
750	50	Yes	-	Cross Valve	Yes	Yes	1/10/80	27510918
750	10	-	-	Cross Valve	Yes	Yes	1/10/80	27510118
750	30	-	Yes	Cross Valve	Yes	Yes	1/10/80	27510218
750	9	-	Yes	Cross Valve	Yes	Yes	1/10/80	27510318
750	9	-	Yes	Lever tap	Yes	Yes	1/10/80	27510418
750	30	-	Yes	Lever tap	Yes	Yes	1/10/80	27510518

* Includes 2 leg straps per shelf box (10 leg bags)

Sterilization
Sterile, Sterilized by Ethylene oxide, Shelf life: 5 years

Certifications	
CE-marking	CE 2797
UKCA-marking	UKCA 0086
Product classification	Class Is (EU)
Quality system certified to	ISO 13485, Medical Device Single Audit Program (MDSAP)
European Medical Device Nomenclature (EMDN)	A06030301 - URINE COLLECTION BAGS
Global Medical Device Nomenclature (GMDN)	58923



Home Care Leg Bags

Sterile Urine Drainage Bags

Instructions for use	
Included?	Yes

Material specifications	
Contains latex?	No
Contains phthalates?	No
Contains cotton?	No
Contains bisphenols?	No
Contains paper?	Packaging only
Bag (including non-return valve)	Polyvinyl chloride
Tube, Connector	Polyvinyl chloride
Cap	Polyethylene
Bottom Outlet	Polyvinyl chloride, Polyethylene, Polypropylene
Tube on Bottom Outlet	Silicone
Straps	Spandex, Nylon, Dacron, Silicone

Packaging information									
Reference	Primary pack./Single peelpack		Secondary pack./Shelf box		Tertiary pack./Shipper			Gross weight kg	Pallet Qty. pcs
	Qty. pcs	Dimensions mm (LxW)	Qty. pcs	Dimensions mm (LxWxD)	Qty. pcs	Dimensions mm (LxWxD)	Volume m ³		
23510118	1	260x190	10	240x190x135	80	560x395x260	0.058	4.50	2,240
23510218	1	260x190	10	240x190x135	80	560x395x260	0.058	4.50	2,240
25010118	1	260x190	10	240x190x135	80	560x395x260	0.058	4.50	2,240
25010218	1	260x190	10	240x190x135	80	560x395x260	0.058	4.50	2,240
25010318	1	260x190	10	240x190x135	80	560x395x260	0.058	4.50	2,240
25010418	1	260x190	10	240x190x135	80	560x395x260	0.058	4.50	2,240
25010518	1	260x190	10	240x190x135	80	560x395x260	0.058	4.50	2,240
25010618	1	260x190	10	240x190x135	80	560x395x260	0.058	4.50	2,240
25010718	1	260x190	10	240x190x135	80	560x395x260	0.058	4.50	2,240
25010818	1	260x190	10	240x190x135	80	560x395x260	0.058	5.51	2,240
27510918	1	260x190	10	240x190x135	80	560x395x260	0.058	5.60	2,240
27510118	1	260x190	10	240x190x135	80	560x395x260	0.058	5.60	2,240
27510218	1	260x190	10	240x190x135	80	560x395x260	0.058	4.50	2,240
27510318	1	260x190	10	240x190x135	80	560x395x260	0.058	4.50	2,240
27510418	1	260x190	10	240x190x135	80	560x395x260	0.058	5.75	2,240
27510518	1	260x190	10	240x190x135	80	560x395x260	0.058	6.40	2,240

Packaging material	
Primary package	Paper, Polyethylene terephthalate/Polyethylene Film
Secondary package	Single-wall corrugated fiberboard
Tertiary package	Double-wall corrugated fiberboard



Home Care Leg Bags

Sterile Urine Drainage Bags

GTIN codes			
REF	Primary package	Secondary package	Tertiary package
23510118	06936881154430	16936881154437	26936881154434
23510218	06936881154447	16936881154444	26936881154441
25010118	06936881154454	16936881154451	26936881154458
25010218	06936881154461	16936881154468	26936881154465
25010318	06936881154478	16936881154475	26936881154472
25010418	06936881154485	16936881154482	26936881154489
25010518	06936881154492	16936881154499	26936881154496
25010618	06936881154508	16936881154505	26936881154502
25010718	06936881154515	16936881154512	26936881154519
25010818	06936881154522	16936881154529	26936881154526
27510918	06936881154539	16936881154536	26936881154533
27510118	06936881154546	16936881154543	26936881154540
27510218	06936881154553	16936881154550	26936881154557
27510318	06936881154560	16936881154567	26936881154564
27510418	06936881154577	16936881154574	26936881154571
27510518	06936881154584	16936881154581	26936881154588

	This way up
	Keep dry
	Fragile, handle with care
	Keep away from sunlight
	Do not re-sterilize
	Do not re-use
	Do not use if package is damaged and consult instructions for use
	Temperature limit (-15° to 49°C)
	Consult instructions for use or consult electronic instructions for use

DISTRIBUTED BY
UNOQUIP GMBH
 Schützengraben 20
 8200 Schaffhausen
 Switzerland
info@unoquip.com
www.unoquip.com

SPONSORED IN AUSTRALIA BY
 UNOQUIP Pty Ltd
 PO Box 1117
 Windsor VIC 3181
 1300 55 80 79
info@unoquip.au



Well Lead Medical Co., Ltd.
 C-4 Jinhu Industrial Estate
 Hualong, Panyu Guangzhou
 511434 P.R. China



Shanghai Intl. Holding
 Corp. GmbH
 Eiffestrasse 80
 20537 Hamburg
 Germany



QS Engineering AG
 Rütihof 1
 CH 8820 Wädenswil
 Switzerland

UK Responsible Person (UKRP)

SUNGO Certification
 Company Limited
 3rd Floor, 70 Gracechurch
 Street, London, EC3V 274HR
 England, United Kingdom

Home Care Leg Bags

Sterile Urine Drainage Bags



Certified Quality Management System
DIN EN ISO 13485:2021
EN ISO 13485:2016 / A11:2021
ISO 13485:2016