

Blue Steam Tape with Indicators

Description

Blue process indicator tape for steam sterilization
 Ref: 31BL1950VE – 31BL2550VE
 Crepe paper coated with a dry natural rubber based adhesive.
 This tape does not include natural rubber latex.
 Tape individually packed with a plastic protective film.

Sterilization

Steam.

Application

Identification and closure for all types of sterilization pack.
 Use a tape dispenser in order to make the cutting easier.

Performances

Good adhesion on wrap made with cellulose fibres.
 Good resistance to sterilization process.
 Indicator colour change : from green to brown/black.

Storage

Products must be stored in their original packaging, in dry and cool areas away from humidity.

Shelf Life

3 years from date of manufacturing

Conformity

ISO 11140-1 (indicator)

Origin

Canada

CE Mark

NA

Technical Datas

PARAMETERS	UNITS	VALUES
<i>THICKNESS</i>	mm	0,17
<i>METAL ADHERENCE</i>	N/25 mm	8
<i>TENSILE STRENGTH</i>	N/25 mm	90
<i>LENGTHENING</i>	%	7
<i>SOLVENT AND EROSION STRENGTH</i>		Good
<i>HEAT AND MOISTURE RESISTANCE</i>		Excellent
<i>HANDLING.</i>		Excellent