




TECHNICAL DATA SHEET

STEAM STERILIZATION

Process indicator Type 1



ADHESIVE LABELS 29X28 MM

PRODUCT CODE: 2928RBT05STVS

PRODUCT DESCRIPTION

Labels with double adhesive support and process indicator for steam sterilization

Label size: 29x28mm

MATERIAL

Label: Vellum paper on double adhesive support

Adhesive: Acrylic

Process indicator: STEAM Standard ISO11140-1 TYPE 1 non toxic ink

Latex free product

PROCESS INDICATOR

STEAM Standard ISO11140-1 Type 1

The process indicator , after the sterilization cycle, turns from pink to brown.

PACKAGING

Each Box contains 30 rolls of 500 labels

Each Carton contains 10 boxes

SHELF LIFE

3 years from the production date

CONDITION OF STORAGE

Keep the product intact into the original package.

Temperature : min. 5°C –max 30°C. Humidity: 35% - max 50%

PRODUCT CONFORMITY

EN ISO 11140-1:2014 Type 1 Indicator

Manufactured in Factory with environmental management system certified EN ISO 14001:2015 and Quality Certification EN ISO 9001:2015