

Masstabelle mediven® thrombexin® 18

Maßtabelle							
Größe	Beinmaße* (Umfang in cm)				Länge	Farbmarkierung	
	Fessel cB	Wade cC	unter Knie cD	Oberschenkel cG		Ferse	Oberschenkel
Knie (AD)					Länge AD		
XS	17 – 19	26 – 32	25 – 31	–	35 – 46	orange	–
S	20 – 22	30 – 38	29 – 37	–	35 – 46	grün	–
M	23 – 25	34 – 44	33 – 43	–	35 – 46	violett	–
L	26 – 28	38 – 48	37 – 47	–	35 – 46	blau	–
L-X	26 – 28	44 – 55	43 – 54	–	35 – 46	grau	–
XL	29 – 31	44 – 55	43 – 54	–	35 – 46	rot	–
XXL	32 – 34	51 – 66	50 – 65	–	35 – 46	braun	–
Oberschenkel (AG)					Länge AG		
XS	17 – 19	26 – 36	25 – 35	40 – 50	60 – 84	orange	orange
S	20 – 22	30 – 42	29 – 41	48 – 60	65 – 89	grün	grün
S-X	20 – 22	30 – 42	29 – 41	56 – 70	65 – 89	grün	grau
M	23 – 25	34 – 46	33 – 45	56 – 70	65 – 89	violett	violett
M-X	23 – 25	34 – 46	33 – 45	64 – 80	65 – 89	violett	grau
L	26 – 28	38 – 51	37 – 50	64 – 80	65 – 89	blau	blau
L-X	26 – 28	38 – 51	37 – 50	72 – 90	65 – 89	blau	grau
XL	29 – 31	42 – 55	41 – 54	72 – 90	65 – 89	rot	rot
XL-X	29 – 31	42 – 55	41 – 54	80 – 100	65 – 89	rot	grau
Hüfte (AG)					Länge AG		
S	20 – 22	30 – 42	29 – 41	48 – 60	65 – 89	grün	grün
M	23 – 25	34 – 46	33 – 45	56 – 70	65 – 89	violett	violett
L	26 – 28	38 – 51	37 – 50	64 – 80	65 – 89	blau	blau
L-X	26 – 28	38 – 51	37 – 50	72 – 90	65 – 89	blau	grau