

Ryles Tube Products

Product Data Sheet

Description	Ryles Tube, Graduated, X-Ray with closed tip, four side eyes, and colour coded funnel connector and attached luer adaptor.
Intended Use	Ryles Tube is designed to provide access to nasogastric tract for introduction of nutrition and aspiration of intestinal secretion. The main use of a Nasogastric tube is for feeding and for administering drugs and other oral agents. A Nasogastric tube is inserted through the nose with the end of the tube reaching into the stomach.
Duration of Use	Short Term Use: for more than 60 minutes and less than 30 days.
Materials	Tube, Funnel Connector and Luer - PVC (DEHP Free); Radio-Opaque Line (X-Ray).
Product Information	Smooth low friction surface of tube facilitates an easy intubation. Soft, frosted and kink resistance PVC tubing. 105cm tube with radio-opaque line provided to its full length and marked at 50, 60, and 70 cm from the tip for accurate placement of tube. Closed distal end with four lateral eyes to reduce risk of blockage. Color coded funnel end connector for easy identification of size and attached luer adaptor. Tube shore hardness: 75~78A. Available in different sizes from FG 6 to FG 20. Supplied sterile, in individual peel pouches.
Quality System	EN ISO 13485:2016
EU Classification	Class IIa
Sterilised by	Ethylene Oxide
Shelf Life	5 years
Storage & Safety	Toxicity: Manufactured from biocompatible materials. Storage: Store in between 5°C to 35°C, avoid excessive heat, protect from direct sunlight and moisture.


Product Code	Tube Length	CH/Dia	Connector	Side Eyes	Tip	Box Qty	
AERT2306	105cm	6ch ø 2.0mm	Funnel + Luer	4	Closed Tip	50	
AERT2308	105cm	8ch ø 2.7mm	Funnel + Luer	4	Closed Tip	50	
AERT2310	105cm	10ch ø 3.3mm	Funnel + Luer	4	Closed Tip	50	
AERT2312	105cm	12ch ø 4.0mm	Funnel + Luer	4	Closed Tip	50	
AERT2314	105cm	14ch ø 4.7mm	Funnel + Luer	4	Closed Tip	50	
AERT2316	105cm	16ch ø 5.3mm	Funnel + Luer	4	Closed Tip	50	
AERT2318	105cm	18ch ø 6.0mm	Funnel + Luer	4	Closed Tip	40	
AERT2320	105cm	20ch ø 6.7mm	Funnel + Luer	4	Closed Tip	40	

*Please note Pennine does not recommend non ENFit tubes for feeding in markets that have adopted the ENFit System.

© 2022 Pennine Healthcare or one of its subsidiaries.
All rights reserved.

The information given in this sheet is not intended to be exhaustive.

The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.

 Poly Medisure Limited, Plot No. PA-010-018 & PA-010-019,
Mahindra World City (JAIPUR) Ltd. Multi-Product SEZ.
Distt- Jaipur 302037 (Rajasthan) India

EC REP Obelis S.A. Boulevard General Wahis 53,
B-1030, Brussels, Belgium. mail@obelis.net

