




CORE Room Turnaround Kits®

Tourniquet Covers

CODE	DESCRIPTION	SIZE	CASE QTY	CASE DIMENSIONS
 50-TQTC01	Tourniquet Cover, 17x17" (44x44cm) w/ Fenestration, 2x2" (5x5cm)	S	200	10x10x9" (25x25x22cm)
 50-TQTC02	Tourniquet Cover, 17x17" (44x44cm) w/ Fenestration, 3x4" (7x10cm)	M	200	10x10x9" (25x25x22cm)
 50-TQTC03	Tourniquet Cover, 23x23" (58x58cm) w/ Fenestration, 5x6" (12x16cm)	L	200	12x12x9" (31x31x22cm)

Materials

-  Kraton Conformable Fenestration Material

Product Attributes

- Conformable fenestrations

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications


- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/Registrations


- EN ISO 13485:2016

Storage


- General warehouse conditions



50-TQTC01



50-TQTC02



50-TQTC03



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-002 2022
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

THINK.MAKE.PROTECT.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384



Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.emergobyul.com

A GRI COMPANY

