

STARFIT™ DRAINAGE BAG SYSTEM



GASTRIC DRAINAGE WITH CONFIDENCE

Measure gastric residual volumes safely and confidently with a drainage system designed with ENFit® connectors – no adapters, no off-label workarounds, just clinical accuracy.



Hoist
medical

Improving quality of life with assisted feeding



STARFIT™ DRAINAGE BAG SYSTEM

DISCOVER CONFIDENCE IN GASTRIC DRAINAGE

The ENFit® Compliant Drainage Bag System

In modern clinical care, **every detail matters** – especially when it comes to enteral access and drainage.

For years, front line caregivers have faced the daily challenge of improvising gastric drainage solutions with urinary collection bags and mismatched connectors. Makeshift setups not only increase risk but also complicate accurate gastric residual volume (GRV) monitoring – a critical component in managing patients with compromised gastric motility.

The **ENFit® Compliant Enteral Drainage Bag System** is engineered with your patient's safety and your workflow in mind. Built to meet the **ISO 80369-3 standard**, this system completely eliminates the risk of misconnections.

No adapters. No workarounds. Just **secure, reliable drainage you can trust –right out of the package.**

Designed for continuous 24-hour use, our system allows for easy, accurate measurement of gastric residuals whenever clinical judgment demands.

The bag hangs conveniently at the bedside, below stomach level, promoting natural gravity drainage while ensuring patient comfort. The integrated hydrophobic air valve allows air to escape. This may also help to relieve gastric bloating, pain, and discomfort often seen in patients with delayed gastric emptying.

Whether your patients are in the ICU, post-operative recovery, or long-term care, this system provides a clinically sound alternative to off-label drainage methods – giving clinicians the confidence to focus on care, not improvisation.



Why Choose Our ENFit® Compliant Drainage Bag System?



Safe ENFit Connection – ISO 80369-3 Compliant



Protects care givers from exposure to gastric fluids



Accurate Gastric Residual Volume (GRV) Monitoring



Use for up to 24 hours



Hydrophobic Valve for Air Venting – Helps Relieve Discomfort



Available in 4 sizes for all clinical needs – neonates to adults



Large bore tubing to facilitate gastric drainage



Bed side hanger included



A True Enteral solution – No Makeshift Adaptors Needed

No more compromising.

Empower your clinical practice with the Enteral Drainage Bag System – the professional choice for safe, effective gastric drainage.



STARFIT™ DRAINAGE BAG SYSTEM



Hoist medical

Improving quality of life with assisted feeding

Give your patients and carers
the confidence they deserve
with an **Enteral Drainage Bag
System** that is designed for
compliance, **ease of
use, reliability, and superior
patient care.**

EXPERIENCE THE DIFFERENCE!

**CONTACT US TODAY
TO LEARN MORE**



hoist-medical.com

**STARFIT™
DRAINAGE BAG SYSTEM**



Sales and General Enquiries
+33(0) 4 91 76 73 29
commercial@hoist-medical.com