

REF 108767AA_(1303245)_Drape Central Venous
Product specifics

- Dimension is 90x75cm
- 45 x 45cm Ø 12,5cm + ZIP + tape
- Zipper feature; without using a scissor the drape can be easily eliminated

Product reference codes

	Item number (Non-sterile)	REF (Sterile)	REF (Non sterile)
Product reference	1303245	108767AA	On request
UDI	08710685109326	08710685035182	On request
Classification	Class I	Class Is	Class I
Notified body	Not applicable	BSI (2797)	Not applicable
Sterilization	To be sterilized before use	Product is sterilized by ethylene oxide (EO)	To be sterilized before use


Intended purpose

The drapes, gowns and equipment covers are intended to minimize the spread of micro-organisms, herewith reducing the risk for post-operative wound infection.

Regulatory information

Product CE marked according MDR 2017/745

Product is EN13795 compliant (test results available upon request)

Dimensions, material use

● 3-ply ○ Transparent PE

Manufacturer

Medica Europe BV
Galliërsweg 20
5349 AT OSS
The Netherlands
www.medica-europe.nl