

TECHNICAL DATA SHEET



Description

Reference

69600: Impervious Comfort Gown, Knit Cuffs, Large
69601: Impervious Comfort Gown, Knit Cuffs, XX-Large
69602: Impervious Comfort Gown, Thumb Hooks, Large
69603: Impervious Comfort Gown, Thumb Hooks, XX-Large

Impervious Film Laminate- Blue, Disposable Fluid-Resistant, Potential for Soaking and Heavy Fluid Contact in Front, Perforated Open Back.

Dimensions

	L	XXL
Length (from shoulder to hem)	121 cm	121 cm
Sleeve length (underarm to cuff/hook)	61 cm	62 cm
Neck length	27 cm	27 cm
Body width (chest)	71 cm	84 cm
Width (overall) :	119 cm	140 cm

Properties

Basis weight of base material: 38 g/m²
Material was tested as per ASTM F 1671- Viral Penetration and as per ASTM F 739-99a- Resistance to Permeation by Liquids and Gases.

Indication

To prevent the transfer of infective agents between the operating room personnel and the patient.
Single use product.

Counter indication


None in particular.

Main materials

Base material: Spunbond-Film Laminate (polypropylene)
Cuffs and Threads: Polyester
Neck Closure: Fastening Tape with Release Liner.
Does not contain Natural rubber latex.
Does not contain DEHP.

Sterilisation

Products are non-sterile.



Packaging

69600, 69601, 69602, and 69603: Shipping case of 100 units.
10 units are placed within 1 polyethylene bag,
10 polyethylene bags are placed within 1 shipping case.
Bar coding: GS1-128 symbology, linear, on shipping case and bags.

Manufacturing

Products are manufactured in China.
The quality system of the manufacturing sites is ISO 13485 compliant.

Regulatory information

Products CE marked as per 93/42/EEC Directive on Medical Devices.
Class of the device: I.

Storage

Store in a dry and cool place, away from intense sources of heat.
Keep as much as practicably possible in its shipper box.

Shelf life

Non-specified.