

Dust Cover bags

DESCRIPTION

Protective packaging, transparent both sides with a self-seal closure (adhesive).

Article code : DC*****

2 sides made of 50µ polypropylene printed with 1 colour (black).

Self-seal closure (adhesive).

Insert allowing writing information on the content

STERILIZATION

Dust Cover bags must not be sterilized

APPLICATION

Protective packaging designed to prevent damage to the sterile barrier system and/or the packaging system, its content during the transport and the storage (according to EN ISO 11607-1)

PERFORMANCES

- Excellent mechanical strength.
- High quality adhesive allowing to reopen the bag numerous times

STORAGE

It is recommended that the products are kept in the original, closed transport carton and are stored in dry and clean conditions protected from direct sunlight and excessive moisture.

2 years from the manufacturing date, provided the above storage conditions are met.

CONFORMITY

Protective packaging according to EN ISO 11607-1

ORIGIN

France

Technical Data

PROPERTIES	UNIT	TOLERANCE	VALUE
MASS	g/m ²	±5%	40,5/45
MELTING POINT	C°	±3	146
BRILLIANCE (GLOSS)	brilliant point	±10	138
TRANSPARENCY	%	±1	2,5
RUBBING COEFFICIENT	/	±0,05	0,1
TREATED SIGHT	/	±0,05	/
NON TREATED SIGHT	/	±0,05	/
LIMIT	N/mm ²	±3	15
ELASTICITY	N/mm ²	±3	14
BREAKING RESISTANCE SE	N/mm ²	±10	40
BREAKING RESISTANCE ST	N/mm ²	±10	30
TENSILE STRENGTH SE	%	±100	700
TENSILE STRENGTH ST	%	±100	700
PERMEABILITY H2O WHEN 38C°	g/m ² /24h cc/m ² /24h	±1	6
PERMEABILITY O2 WHEN 25C°	Bar	±200	1600
SEALING STRENGTH	N/mm ²	±8	28