

Product Specification for Thermo Mini Poncho MP101


Details

	Reference:
Product standards	EN 13795-1:2019
Classification according to MDR	Class I, Medical Device Regulation (EU) 2017/745
Certification of Manufacturer	ISO 9001:2015 CN14/21396 ISO 13485:2016 and EN ISO 13485: 2016 CN14/21395
Environment Management	ISO 14001 CN23919E10173R0S

Basic UDI-DI

574400129-MP101-2S

Physical Dimensions in cm and Properties



Dimension	A	B	C	D	E	F	G
Specification	20+/- 0,5cm	14,5+/-0,5 cm	45+/-1cm	24+/-1cm	43+/-1cm	17+/-0,5cm	5+/-0,5cm
Units:	Colour:	Material:				Shelf Life:	Average Weight:
Cm	Dark Blue	Polypropylene Spunbond – and Composite non-woven Fabric				3 Years	47 gsm +/- 10%