

UnoMeter™ Safeti Plus

Hourly Diuresis Monitoring System



Product range

UnoMeter™ Safeti Plus – with KombiKon™ sample port. Hourly Diuresis Monitoring System

REF No	Description	Chamber volume	Bag volume	Sample Port	Double lumen tube	NRV in bag	NRV in inlet tube	Bottom Outlet (free flow)	Round Hook	Straps	Cords	Exchangeable bag	Sliding clamp
1710086	2.0L Exchangeable bag, straps and round hook, 110cm tube, Free flow bottom outlet, sterile	500 ml	2000 ml	KombiKon™ Needle-free sample port with inspection window	●	●	●	●	●	●		●	●
1710085	2.0L Exchangeable bag, straps and round hook, 150cm tube, Free flow bottom outlet, sterile	500 ml	2000 ml	KombiKon™ Needle-free sample port with inspection window	●	●	●	●	●	●		●	●
1710088	2.0L Exchangeable bag, straps and cord, 110cm tube, Free flow bottom outlet, sterile	500 ml	2000 ml	KombiKon™ Needle-free sample port with inspection window	●	●	●	●		●	●	●	●
1710087	2.0L Exchangeable bag, straps and cord, 150cm tube, Free flow bottom outlet, sterile	500 ml	2000 ml	KombiKon™ Needle-free sample port with inspection window	●	●	●	●		●	●	●	●
1710089	2.0L Exchangeable bag, straps and round hook, 120cm tube, Free flow bottom outlet, sterile	500 ml	2000 ml	KombiKon™ Needle-free sample port with inspection window	●	●	●	●	●	●		●	●
1710090	2.0L Exchangeable bag, straps and round hook, 200cm tube, Free flow bottom outlet, sterile	500 ml	2000 ml	KombiKon™ Needle-free sample port with inspection window	●	●	●	●	●	●		●	●

Material specifications:

Component	Materials
Contains latex or DEHP?	No
KombiKon™ connector and cap	PVC – Phthalate free, TPE, PP, PE
Tube, Clamp and anti-kink spring	PVC – Phthalate free, PP, PE, TPU
Rigid chamber incl. strap, Hanger, lid and inner cylinder	SAN, ABS, ePTFE membrane, PE, SBC, PA, PVC – Phthalate free, NBR
Bag incl. bottom outlet	PVC – Phthalate free, ePTFE membrane, PE, ABS
Packaging materials	PE/PET, PE, corrugated cardboard

Sterilization:

Sterile?	Yes
Sterilized by	Ethylene oxide
Shelf life	3 years

Certifications:

CE-marking?	Yes
Product classification	Class Is, Im (EU)
Quality system certified to	ISO 13485
Global Medical Device Nomenclature (GMDN)	32072
EMDN Code	A06030302

Storage and handling

Instructions for use included.



Packaging specifications:

REF No	Shipping Carton				
	Quantity (pcs)	Dimensions (mm) L x W x H	Gross Weight (kg)	Piece identifier GTIN - UDI	Inner carton identifier GTIN - UDI
1710085	14	580x380x280	6,8	07350057674253	07350057674260
1710086	14	580x380x280	6,5	07350057674284	07350057674291
1710087	14	580x380x280	6,6	07350057674314	07350057674321
1710088	14	580x380x280	6,3	07350057674345	07350057674352
1710089	14	580x380x280	6,6	07350057674376	07350057674383
1710090	14	580x380x280	7,0	07350057674406	07350057674413

