

# Uterine Aspiration Curette

## Product Data Sheet

<b>Description</b>	Uterine Aspiration Curette
<b>Intended Use</b>	Uterine Suction Catheter connected to uterine suction tube is used for curettage and evacuation of the endometrial contents during induced termination of pregnancy or after miscarriage or as a curative measure in cases of dysmenorrhea.
<b>Duration of Use</b>	Normally intended for continuous use of less than 60 minutes, but will be patient specific depending on the duration of the surgical procedure.
<b>Materials</b>	Curette - K-Resin Tube – PE Connector - ABS
<b>Product Information</b>	The curette comes in sizes ranging from 6 mm to 16 mm in diameter and size 4 mm and 5 mm are additionally available in Flexible Curettes, offering a convenient and adaptable system for the extraction of the contents of the uterus; Curved and straight design for options in Rigid Curettes; Moulded markings for sizing and correct tip orientation in Rigid Curettes; Printed markings for sizing and correct tip orientation in Flexible Curettes; Leak proof seals; Easily connect to rotary swivel handle and wide bore, high vacuum tubing; Sterile, double wrapped, individual peel pouches
<b>Quality System</b>	ISO 13485
<b>EU Classification</b>	Class IIa
<b>Sterilised by</b>	Ethylene Oxide
<b>Shelf Life</b>	5 years
<b>Storage &amp; Safety</b>	Toxicity: Manufactured from biocompatible materials Storage: Store in a cool and dry place, avoid direct sunlight



© 2023 Pennine Healthcare or one of its subsidiaries.  
All rights reserved.

The information given in this sheet is not intended to be exhaustive.

The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.



GCMEDICA ENTERPRISE LTD.(WUXI) Loujin Industrial Park,  
Shuofang, Wuxi 214143, Jiangsu, P.R. China



Luxus Lebenswelt GmbH  
Kochstr.1, 47877, Willich, Germany



# Uterine Aspiration Curette

## Product Data Sheet

## SURGICAL SUCTION

Product Code	Type	Size	Box Qty
AEKC1506	Straight	6mm	50
AEKC1507	Straight	7mm	50
AEKC1508	Straight	8mm	50
AEKC1509	Straight	9mm	50
AEKC1510	Straight	10mm	50
AEKC1511	Straight	11mm	50
AEKC1512	Straight	12mm	50
AEKC1006	Curved	6mm	50
AEKC1007	Curved	7mm	50
AEKC1008	Curved	8mm	50
AEKC1009	Curved	9mm	50
AEKC1010	Curved	10mm	50
AEKC1011	Curved	11mm	50
AEKC1012	Curved	12mm	50
AEKC1514	Straight, with connector	14mm	30
AEKC1014	Curved, with connector	14mm	30
AEAC3004	Flexible	4mm	50
AEAC3005	Flexible	5mm	50
AEAC3006	Flexible	6mm	50
AEAC3007	Flexible	7mm	50
AEAC3008	Flexible	8mm	50
AEAC3009	Flexible	9mm	50
AEAC3010	Flexible	10mm	50
AEAC3011	Flexible	11mm	50
AEAC3012	Flexible	12mm	50

© 2023 Pennine Healthcare or one of its subsidiaries.  
All rights reserved.

The information given in this sheet is not intended to be exhaustive.

The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.



GCMEDICA ENTERPRISE LTD.(WUXI) Loujin Industrial Park,  
Shuofang, Wuxi 214143, Jiangsu, P.R. China



Luxus Lebenswelt GmbH  
Kochstr.1, 47877, Willich, Germany

