

Veress Needles

Disposable Veress Needles



Veress Needle




K-Line Veress Needle

Veress Needles


The Veress Needle is a sterile single patient use product to create a pneumoperitoneum prior to laparoscopic procedures.

- Sharp bevelled needle tip ensures simple abdominal wall penetration
- Transparent handle with clear red visual position indicator
- Spring-loaded stylet provides audible and tactile feedback of correct needle placement
- Blunt stylet shields needle tip once correct positioning has been achieved
- High-flow rotating stopcock with luer lock connector for the connection with the insufflation tube
- Available in two lengths: 120mm and 150mm
- Available both as sterile single-packed product and for inclusion in Mölnlycke® Procedure trays

Veress Needle




Veress Needle body




Veress Needle tip

K-Line Veress Needle



Veress Needle body



Veress Needle tip

Ref. No.	Product name and description	Size mm	Pcs/Disp
899105	Veress Needle	120	25
899106	Veress Needle	150	25
800105	K-Line Veress Needle	120	10
800106	K-Line Veress Needle	150	10

Legal manufacturer: Unimax Medical Systems Inc.
Legal manufacturer: Hangzhou Kangji Medical Instrument Co., Ltd.
Imported by: Mölnlycke Health Care AB

Find out more at [molnlycke.com](https://www.molnlycke.com)

Mölnlycke Health Care AB, P.O. Box 6, Entreprenörstråket 21, SE-431 21 Mölndal Sweden. Phone + 46 31 722 30 00.
The Mölnlycke and ProcedurePak trademarks, names and logos are registered globally to one or more of the Mölnlycke Health Care group of companies. ©2025 Mölnlycke Health Care AB. All rights reserved. HQIM007666

